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  BIO-LA-TEST®  																				
CREATININE (CREAT 100)																				
Cat. No. 10003145 Store at (+2 to +8)°C																				
Reagent Set for the determination of creatinine in serum or in urine. The Set is sufficient for at least 100 analyses.																				
Principle Creatinine reacts in an alkaline medium with picric acid to form an orange-red coloration, which is suitable for photometric determination. Serum creatinine is determined after deproteinizing the sample, the urine is analyzed after dilution with water. The method is unspecific and compounds with active methylene group and some of the reducing substances (glucose) interfere with the determination.																				
References Löken, F.: Scand. J. Clin. Lab. Invest. 12, 381 (1954) van Pilsen, J. F., Boris, M.: Clin. Chem. 3, 90 (1957)																				
Reagents <table style="width: 100%; border: none;"> <tr> <td style="width: 80%;">1 Creatinine standard</td> <td style="text-align: right;">(1 vial)</td> </tr> <tr> <td>creatinine 442.5 µmol/l</td> <td></td> </tr> <tr> <td>2 Albumin</td> <td style="text-align: right;">(1 vial)</td> </tr> <tr> <td>albumin 0.16 g/vial</td> <td></td> </tr> <tr> <td>3 Deproteinizing solution</td> <td style="text-align: right;">(55 ml)</td> </tr> <tr> <td>trichloroacetic acid 1.22 mol/l</td> <td></td> </tr> <tr> <td>4 Picric acid</td> <td style="text-align: right;">(55 ml)</td> </tr> <tr> <td>solution of 0.04 mol/l</td> <td></td> </tr> <tr> <td>5 Sodium hydroxide</td> <td style="text-align: right;">(55 ml)</td> </tr> <tr> <td>solution of 0.75 mol/l</td> <td></td> </tr> </table>	1 Creatinine standard	(1 vial)	creatinine 442.5 µmol/l		2 Albumin	(1 vial)	albumin 0.16 g/vial		3 Deproteinizing solution	(55 ml)	trichloroacetic acid 1.22 mol/l		4 Picric acid	(55 ml)	solution of 0.04 mol/l		5 Sodium hydroxide	(55 ml)	solution of 0.75 mol/l	
1 Creatinine standard	(1 vial)																			
creatinine 442.5 µmol/l																				
2 Albumin	(1 vial)																			
albumin 0.16 g/vial																				
3 Deproteinizing solution	(55 ml)																			
trichloroacetic acid 1.22 mol/l																				
4 Picric acid	(55 ml)																			
solution of 0.04 mol/l																				
5 Sodium hydroxide	(55 ml)																			
solution of 0.75 mol/l																				
Composition of reaction mixture <table style="width: 100%; border: none;"> <tr> <td style="width: 80%;">Trichloroacetic acid</td> <td style="text-align: right;">0.15 mol/l</td> </tr> <tr> <td>Picric acid</td> <td style="text-align: right;">10.00 mmol/l</td> </tr> <tr> <td>Sodium hydroxide</td> <td style="text-align: right;">0.19 mol/l</td> </tr> <tr> <td>Volume ratio of serum/reaction mixture</td> <td style="text-align: right;">1/8</td> </tr> </table>	Trichloroacetic acid	0.15 mol/l	Picric acid	10.00 mmol/l	Sodium hydroxide	0.19 mol/l	Volume ratio of serum/reaction mixture	1/8												
Trichloroacetic acid	0.15 mol/l																			
Picric acid	10.00 mmol/l																			
Sodium hydroxide	0.19 mol/l																			
Volume ratio of serum/reaction mixture	1/8																			
Reference values fS Creatinine (µmol/l) Men 61 – 115 Women 53 – 97 The range of reference values is only approximate, it is recommended to all laboratories to verify the extension of reference interval for their concrete examined population.																				
Calibration and quality control BIO-LA-TEST® CREATININE STANDARD, Cat. No. 10003146																				

Working solutions

- Solution 1** Dissolve the content of the vial with Reagent 1 in 8.00 ml of distilled water. Solution contains creatinine 442.5 µmol/l.
 Stability: at least 1 week when stored at (+2 to +8)°C.
- Solution 2** Dissolve the content of the vial with Reagent 2 in 8.00 ml of distilled water. Solution contains albumin 20 g/l.
 Stability: at least 1 week when stored at (+2 to +8)°C.
- Solution 3** By mixing 4.00 ml of Solution 1 with 6.00 ml of Solution 2 prepare standard solution with the content of creatinine 177 µmol/l.
 Stability: at least 1 week when stored at (+2 to +8)°C.

Procedure

Samples: non-haemolytic serum, heparinized or EDTA plasma, urine

Wavelength (500–510) nm
 Cuvette 1 cm
 Temperature (+15 to +25)°C

Determination in serum

Pipette (ml)	Sample	Control solution	Standard
Serum	0.50	–	–
Dist. water	1.00	1.50	1.00
Solution 3	–	–	0.50
Reagent 3	0.50	0.50	0.50
Stir, centrifuge and into another test tubes pipette.			
Supernatant	1.00	1.00	1.00
Reagent 4	0.50	0.50	0.50
Reagent 5	0.50	0.50	0.50
Stir and after exactly 20 min read the absorbances of sample (A ₁) and standard (A ₂) against control solution.			

Determination in urine

Dilute the urine sample with distilled water in the ratio 1+100 prior to use (result x 50).

Pipette (ml)	Sample	Control solution	Standard
Diluted Urine	0.50	–	–
Solution 3	–	–	0.25
Dist. water	0.25	0.75	0.50
Reagent 3	0.25	0.25	0.25
Reagent 4	0.50	0.50	0.50
Reagent 5	0.50	0.50	0.50
Stir and exactly after 20 min read the absorbances of sample (A ₁) and standard (A ₂) against control solution.			

Calculation

$$\text{Creatinine } (\mu\text{mol/l}) = 177 \cdot \frac{A_1}{A_2}$$

Notes

If the content of creatinine in the sample exceeds 440 µmol/l, dilute the sample with distilled water and repeat the analysis (result x dilution).

Sodium hydroxide absorbs carbon dioxide. Keep the bottle with the Reagent 5 closed tightly.

Performance characteristics

Measured on Hitachi 911, without deproteinisation, calibrated on value „Blank rated and compensated“ (compensation of bilirubin and haemoglobin interferences).

Linearity: up to 1 mmol/l, Detection limit: 4.59 µmol/l,

Measuring range: 13.95 - 1000 µmol/l

Precision (at 37°C)

INTRA-ASSAY n = 20	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	123.19	2.37	1.92
Sample 2	370.58	3.43	0.93
Sample 3	1146.06	8.36	0.73
INTER-ASSAY n = 10	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	122.93	2.93	2.38
Sample 2	368.98	5.84	1.58
Sample 3	1144.86	11.96	1.04

Accuracy

Test was compared with commercially available test. Results:

N = 100, r = 1.000, y = 0.932 x + 11.86 µmol/l

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person. Reagent 3 contains trichloroacetic acid (20%), which is corrosive and dangerous for environment. Reagent 5 contains corrosive sodium hydroxide (3.7%).

Reagent 3

R35 Causes severe burns.
 R 51/53 Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

S 24 Avoid contact with skin.

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S 36/37/39 Wear suitable protective clothing, gloves and eye/face protection.

S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

S 60 This material and its container must be disposed of as hazardous waste.

S 61 Avoid release to the environment. Refer to special instructions/Safety data sheets.

Reagent 5

R34 Causes burns.

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S37/39 Wear suitable gloves and eye/face protection.

S45 In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

While working it is especially significant to observe sanitary rules, not to eat, drink, or smoke and to use self preventatives.

First aid: At an accidental ingestion drink up approx. 0.5 l of water. At an eye contact flush immediately and thoroughly with large quantity of water. At a skin contact wash the affected spot with soap and warm water. In all serious cases of health damage consult the physician.

Waste disposal: All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging paper waste and rinsed containers to recycling.

Date of last revision: 27. 2. 2013

Produced	Erba Lachema s.r.o. Karásek 1d, 62100 Brno, CZ e-mail: diagnostics@erbalachema.com , www.erbalachema.com
-----------------	---