BIO-LA-TEST®

**GLUCOSE GOD 250 ES****(GLU 250 ES)**

Cat. No 10003257

Store at  
(+2 to +8)°C

For 250 ml of GOD-POD reagent for enzymatic determination of glucose.

**Principle**

Glucose oxidase catalyses oxidation of glucose to hydrogen peroxide and gluconate. In the presence of peroxidase, hydrogen peroxidase reacts with the chromogen and 4-aminophenazone to form coloured product.

**References**

Trinder, P.: Ann. Clin. biochem. 6, 24 (1969)  
 Chromý, V., Breinek, P., Roháček, J.: Biochem. Clin. Bohemoslov. 16, 275 (1987)  
 Roháček, J.: Chromý, V., Breinek, P.: Biochem. Clin. Bohemoslov. 16, 183 (1987)

**Reagents**

- |   |          |
|---|----------|
| 1 Enzymes                                       | (1 vial) |
| POD ≥4.17 µkat, GOD ≥50 µkat                    |          |
| 4-aminoantipyrine – 0.245 mmol/vial             |          |
| 2 Buffer-Chromogen                              | (250 ml) |
| phosphates 0.14 mol/l, 3-methylphenol 10 mmol/l |          |
| 3 Glucose standard                              | (20ml)   |
| solution 10 mmol/l                              |          |

**Incubation mixture**

Phosphate bufer, pH 8 (25°C)	140 mmol/l
GOD	≥ 166 µkat/l
POD	≥ 16 µkat/l
3-Methylphenol	10 mmol/l
4-Aminoantipyrine	1 mmol/l
Volume ratio serum/incubation mixture	1/101

**Reference values**

cB Glucose (mmol/l)	3.9–5.6
fS, fP Glucose (mmol/l)	3.88–5.83
dU Glucose (mmol/24 hrs)	1.4

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

**Calibration and quality control**

For calibration we recommend  
 BIO-LA-TEST® LYONORM CALIBRATOR, Cat. No. 10003200  
 For control we recommend:  
 BIO-LA-TEST® LYONORM HUM N, Cat. No. 10003204  
 BIO-LA-TEST® LYONORM HUM P, Cat. No. 10003206

**Additional reagents** (not included in the test)

Trichloroacetic acid, solution 50 g/l.

**Working solution**

Dissolve the contents of vial 1 in approx. 6 ml of dist. water and add it into the bottle 2. Mix well.

Stability: 1 month at (+2 to +8)°C in dark,  
 10 days at (+15 to 25)°C in dark.

**Procedure**

Samples: non-haemolytic serum, heparinized or EDTA plasma  
 Wavelength (480–540) nm, maximum at 498 nm  
 Cuvette 1 cm  
 Temperature (+15 to +25)°C or (37 ±0.1)°C

**Assay without deproteinization**

Mix in three test tubes working solution in the ratio 100+1 or 150+1 with serum or plasma (sample), with standard solution of glucose

(standard) and with dist. water (blank). Incubate for 15 min at 37°C or 30 min at (+15 to +25)°C. Avoid exposure to direct light. Read the absorbances of sample (A<sub>1</sub>) and standard (A<sub>2</sub>) against the blank within 30 min after the incubation.

**Assay with deproteinization**

Deproteinize the blood, non-coagulating blood, haemolytic serum or plasma with the solution of trichloroacetic acid. Pipette into centrifuge tube 0.50 ml of deproteinizing solution and 0.05 ml of a sample, shake and centrifuge. Dilute the standard solution of glucose with dist. water in the same way.

Mix in three test tubes working solution in the ratio 10+1 with the supernatant (sample), standard glucose solution diluted with water (standard) and dist. water (blank). Proceed as described under the assay without deproteinization.

**Analyse of urine**

Dilute urine with dist. water 1+10 and proceed as described under the assay without deproteinization (result x 11).

**Calculation**

$$\text{Glucose (mmol/l)} = \underline{a} \cdot \frac{A_1}{A_2}$$

where  $\underline{a}$  is the concentration of glucose in the standard solution in mmol/l.

**Notes**

Vial with Reagent 2 and working solution may contain phosphate crystals which can occur at the low temperature. Redissolve the crystals by gentle heating to approx. 37°C.

If the glucose concentration exceeds 22 mmol/l, repeat the analysis with the sample diluted 1+2 or 1+10 with dist. water (result x dilution).

Anticoagulating agents K<sub>2</sub>EDTA, citrate, oxalate and heparin may be used. The blood should be deproteinized before the determination of glucose. CAUTION! For diluting glucose standard and blank use only dist. water!

Measurement after only 10 min at 37°C is possible in single assay provided that sample and standard are prepared and measured simultaneously.

For kinetic procedure, read the absorbances of the sample and standard (ΔA/min, 37°C) between 30th and 90th seconds.

**Recalculation of units**

1 kat/l = 1 mol/l.s  
 1 µkat/l = 60 U/l

**Stability and storage**

Store reagent at (+2 to +8)°C before and after opening in dark.

Expiration time is mentioned on the cover.

After using the vial close immediately to prevent the evaporation or contamination of the reagent.

**Health protection**

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Any dangerous reagents are not included in kits. If swallowed it is necessary to drink of about 0.5 litre of water. On skin or eye contact rinse the place with jet of tap water.

If necessary, consult a physician.

**Waste disposal**

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and rinsed containers to recycling.

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