CREATININE ENZYMATIC

Cat. No.	Pack Name	Packaging (Content)
BLT00065	CREA ENZ 204	R1: 3 x 50 ml, R2: 3 x 18 ml

EN

CE IVD

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Creatinine in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

Creatinine is a waste product formed in muscle from the high energy storage compound, creatine phosphate. The amount of creatinine produced is fairly constant (unlike Urea) and is primarily a function of muscle mass. It is not greatly affected by diet, age, sex or exercise. Creatinine is removed from plasma by glomerular filtration and then excreted in urine without any appreciable resorption by the tubules.

Creatinine is used to assess renal function, however, serum creatinine levels do not start to rise until renal function has decreased by at least 50%.

PRINCIPLE

In the first reaction, creatinase and sarcosine oxidase are used in the enzymatic hydrolysis of endogenous creatine to produce hydrogen peroxide, that is eliminated by catalase.

Creatinase and 4-aminoantipyrine are added, and only the creatine generated from creatinine by creatinase is hydrolysed sequentially by creatinase and sarcosine oxidase to produce hydrogen peroxide. This newly formed hydrogen peroxide is measured in a coupled reaction catalysed by peroxidase, with N-ethyl-N-sulphopropyl-m-toluidine (ESPMT) as a chromogen.

The absorbance of the produced complex at 546 nm is proportional to the creatinine concentration in the sample.

REAGENT COMPOSITION

R1

Good's buffer pH 7.5	25 mmol/l
Creatinase	12 kU/l
Sarcosine oxidase	8 kU/l
Ascorbate oxidase	2 kU/l
Catalase	200 kU/l
ESPMT	0.47 mmol/l
Detergent	< 1 %
Gentamicin	< 0.1 %

R2

Good's buffer pH 7.5	100 mmol/l
Creatinase	300 kU/l
Peroxidase	20 kU/l
4-aminoantipyrine	2.95 mmol/l
Detergent	< 0.5 %
Sodium azide	< 0.1 %

REACTION MIXTURE

Good's buffer pH 7.5	43.03 mmol/l
Creatinase	8.85 kU/l
Sarcosine oxidase	5.90 kU/l
Ascorbate oxidase	1.48 kU/l
Catalase	47.5 kU/l
ESPMT	0.35 mmol/l
Creatinase	73.77 kU/l
Peroxidase	4.92 kU/l
4-aminoantipyrine	0.73 mmol/l

REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

After first opening, reagents are stable for 30 days at 2–8°C if stored at appropriate conditions, closed carefully and without any contamination.

SPECIMEN COLLECTION & HANDLING

Use serum, plasma (heparin, EDTA), urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum / plasma:	7 days	at 4–25°C
	at least 3 months	at -20°C
in urine:	2 days	at 20–25°C
	6 days	at 4–8°C
	6 months	at -20°C

For the determination in urine use 24 hours specimen. It is important to exactly measure the volume of collected urine. Dilute urine samples in 1+19 ratio with distilled water and multiply results by 20.

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator Lyonorm Calibrator is recommended.

QUALITY CONTROL

For quality control, it is recommended to use following materials:

Lyonorm HUM N, Lyonorm HUM P.

UNIT CONVERSION

mg/dl x 88.4 = µmol/l

EXPECTED VALUES

fS Creatinine (µmol/l)	
male	44 – 97
female	44 – 80
dU Creatinine (mmol/24 hrs)	5 – 18

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 3.7 µmol/l

Linearity: 5700 µmol/l

Measuring range: 3.7 – 5700 µmol/l

Intra-assay precision Within run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	89.83	1.067	1.19
Sample 2	326.89	3.087	0.94

Inter-assay precision Run to run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	89.82	1.99	2.18
Sample 2	313.8	5.40	1.72

COMPARISON

A comparison between XL-Systems CREA ENZ (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 1.004 x + 0.300 \mu\text{mol/l} \quad r = 0.996$$

INTERFERENCES

Following substances do not interfere:

bilirubin up to 30 mg/dl, triglycerides up to 1000 mg/dl, haemoglobin up to 5 g/l.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagents of the kit are not classified as dangerous but Reagent R2 contains less than 0.1% sodium azide – classified as toxic and dangerous substance for the environment.

WASTE MANAGEMENT

Please refer to local legal requirements.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

PROCEDURE

Wavelength: 546 (530–560) nm

Cuvette: 1 cm

Temperature: 37 °C

Serum/reaction mixture ratio 1/61

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

Two-reagent method

	Reagent blank	Calibrator (Standard)	Sample
Reagent R1	0.90 ml	0.90 ml	0.90 ml
Sample	-	-	0.02 ml
Calibrator (Standard)	-	0.02 ml	-
Distilled water	0.02 ml	-	-

Mix and incubate 3-5 min. at 37 °C. Measure absorbance A1 of the sample (A_{sam}), calibrator (A_{cal}) and blank (A_{bl}). Then add:

Reagent R2	0.30 ml	0.30 ml	0.30 ml
------------	---------	---------	---------

Mix and incubate 5-10 min. at 37 °C. Measure absorbance A2 of the sample (A_{sam}), calibrator (A_{cal}) and blank (A_{bl}).

Calculate absorbance change $\Delta A = A2 - A1$ for sample, calibrator and blank.

CALCULATION

$$\text{Creatinine } (\mu\text{mol/l}) = \frac{\Delta A_{sam} - \Delta A_{bl}}{\Delta A_{cal} - \Delta A_{bl}} \times C_{cal}$$

C_{cal} = calibrator (standard) concentration

Applications for automatic analysers will be supplied on request.

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