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SOLUNORM

Kits SOLUNORM are standard stabilized solutions for daily calibration and quality control in clinical chemistry for the analysis of biological fluids. The kits SOLUNORM always contain two types of standard solutions in 1 ml vials. Their concentrations enable calibration and quality control in the range of normal and pathological values. One kit contains two sets of solutions, each of 10x1 ml.

Use the standard solution immediately after opening a vial. Do not store the solution already opened.

Survey of the SOLUNORM kits

Type of the kit	Abbreviation	Solution No.	Contents of 10x1 ml	Cat. No.
SOLUNORM	SONO G	1	glucose: 5 mmol/l ±0.1 mmol/l	10003179
Glucose		2	glucose: 10 mmol ±0.2 mmol/l/l	
SOLUNORM	SONO G U	1	glucose: 5 mmol/l ±0.1 mmol/l	10003178
Glucose – Urea		2	urea: 5 mmol/l ±0.1 mmol/l glucose: 10 mmol/l ±0.2 mmol/l urea: 15 mmol/l ±0.3 mmol/l	
SOLUNORM	SONO Cl	1	sodium chloride: 90 mmol/l ±2 mmol/l	10003173
Chlorides		2	sodium chloride: 110 mmol/l ±3 mmol/l	
SOLUNORM	SONO Na K	1	sodium: 135 mmol/l ±1.3 mmol/l	10003180
Sodium – Potassium		2	potassium: 4 mmol/l ±0.04 mmol/l sodium: 145 mmol/l ±1.3 mmol/l potassium: 5.1 mmol/l ±0.05 mmol/l	
SOLUNORM	SONO OSMO	1	300 mosmol/kg ±3 mosmol/kg	10003182
Osmolality		2	400 mosmol/kg ±4 mosmol/kg	

Kits Solunorm are stored at the temperature of (+2 to +25) °C.

Health protection

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any dangerous reagents are not included in kits. If swallowed it is necessary to drink of about 0.5 litre of water. On skin or eye contact rinse the place with jet of tap water. If necessary, consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging paper waste and containers to recycling.

Date of last revision: 17. 7. 2015

SYMBOLS USED ON LABEL

REF	Catalogue Number	Manufacturer	See Instruction for Use	LOT	Lot Number	CE Mark - Device comply with the Directive 98/79/EC
Storage Temperature	Expiry date	IVD	In vitro Diagnostics	CONT	Content	

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