IMULATEST

MAU STANDARD

Kat. č. 10003230

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Cat. No. 10003230

Store at (+2 to +8) °C

Set of standards for immunoturbidimetric determination of albumin in human urine (3 x 1 ml)

batch	xxxxxx
expiration	mm/rrrr

Certified values

parameter	concentration	uncertainty	unit
Albumin	000	00.0	mg/l

Referential material, where the values are related to:
ERM - DA470k / IFCC

Characteristics

MAU standard is a synthetic liquid calibration material contains human albumin. It is used for calibration during determination of albumin in human urine.

Way of use and stability

The preparation is ready for use. It is treated as standard according to the instruction included to the diagnostic set. If MAU standard is stored correctly at (+2 to +8)°C, it is stable till the expiration date stated on the cover. After use the vials must be screwed tightly and stored at (+2 to +8) °C. Avoid contamination, pipette and decantation of standard back to the vial must be excluded.

Precautions in use

The set is specified only for use in vitro. To be handled by entitled and professionally educated person. Serum was tested for presence of HBsAg and antibodies against HIV and HCV with negative results. In spite of the fact, the presence of some pathogenic materials can not be excluded. For this reason the content of the set has to be treated as potentially infectious material.

Reagents of the kit are not classified like dangerous, but contain in low concentration (<0.1%) sodium azide, very toxic substance, dangerous for the environment.

First aid

If swallowed, a half liter of water should be drunk. In case of contact with skin or eyes, rinse them out with drinking water.
The qualified medical aid has to be provided.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and containers to recycling.

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