IMULATEST

CRP STANDARD (CRPS 3x1)

Kat. č. 10003237

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Cat. No. 10003237

Store at (+2 to +8) °C

Set of standards for immunoturbidimetric determination of C-reactive protein in human serum and plasma (3x 1 ml)

batch	xxxxxx
expiration	mm/yyyy

Certified values

Parameter	Concentration	Uncertainty	Unit
C-reactive protein	000	0.0	mg/l

Referential material, where the values are related to:
ERM-DA474/IFCC

Characteristics

CRP standard is a synthetic liquid calibration material contains human CRP. It is used for calibration during determination of C-reactive protein in human serum and plasma.

Way of use and stability

The preparation is ready for use. It is treated as standard according to the instruction included to the diagnostic set. If CRP standard is stored correctly at (+2 to +8)°C, it is stable till the expiration date stated on the cover. After use the vials must be screwed tightly and stored at (+2 to +8) °C. Avoid contamination, pipette and decantation of standard back to the vial must be excluded.

Precautions in use

The set is specified only for use in vitro. To be handled by entitled and professionally educated person. Serum was tested for presence of HBsAg and antibodies against HIV and HCV with negative results. In spite of the fact, the presence of some pathogenic materials can not be excluded. For this reason the content of the set has to be treated as potentially infectious material.

Reagents of the kit are not classified like dangerous, but contain in low concentration (<0.1%) sodium azide, very toxic substance, dangerous for the environment.

First aid

If swallowed, a half liter of water should be drunk. In case of contact with skin or eyes, rinse them out with drinking water. The qualified medical aid has to be provided.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and containers to recycling.

References

Johnson, A. M.: A new international preparation for proteins in human serum, Arch. Pathol. Lab. Med. 117 (1993) 29-31

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