IMULATEST**Multistandard SET 3****(MSS 3x 1)**

Kat. č. 10003238

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Cat. No. 10003238

Store at (+2 to +8) °C

Set of standards for immunoturbidimetric determination in human serum and plasma (3 x 1 ml)

batch expiration	xxxxxx mm/yyyy
------------------	-------------------

Certified values

Parameter	Concentration 1	Uncertainty	Unit
Albumin	00.00	0.000	g/l
α_1 -Antitrypsin	0.00	0.00	g/l
α_1 acid-glycoprotein	0.00	0.00	g/l
C3 Complement	0.00	0.00	g/l
C4 Complement	0.000	0.000	g/l
Ceruloplasmin	0.000	0.000	g/l
Haptoglobin	0.00	0.00	g/l
Imunoglobulin A	0.00	0.00	g/l
Imunoglobulin G	0.00	0.000	g/l
Imunoglobulin M	0.00	0.00	g/l
Transferrin	0.00	0.00	g/l
Parameter	Concentration 2	Uncertainty	Unit
Albumin	00.00	0.00	g/l
α_1 -Antitrypsin	0.00	0.00	g/l
α_1 acid-glycoprotein	0.00	0.00	g/l
C3 Complement	0.00	0.00	g/l
C4 Complement	0.000	0.000	g/l
Ceruloplasmin	0.000	0.000	g/l
Haptoglobin	0.00	0.00	g/l
Imunoglobulin A	0.00	0.00	g/l
Imunoglobulin G	00.00	0.00	g/l
Imunoglobulin M	0.00	0.00	g/l
Transferrin	0.00	0.00	g/l
Parameter	Concentration 3	Uncertainty	Unit
Albumin	000.00	0.00	g/l
α_1 -Antitrypsin	0.00	0.00	g/l
α_1 acid-glycoprotein	0.00	0.00	g/l
C3 Complement	0.00	0.00	g/l
C4 Complement	0.00	0.00	g/l
Ceruloplasmin	0.00	0.00	g/l
Haptoglobin	0.00	0.00	g/l
Imunoglobulin A	0.00	0.00	g/l
Imunoglobulin G	00.00	0.00	g/l
Imunoglobulin M	0.00	0.00	g/l
Transferrin	0.00	0.00	g/l

Referential material, where the values are related to:
ERM-DA470k/IFCC**Characteristics**

Multistandard set 3 is a liquid calibration material made of human serum mixture. It is used for calibration during determination of common proteins in human serum and plasma.

Way of use and stability

The preparation is ready for use. It is treated as standard according to the instruction included to the diagnostic set. If Multistandard set 3 is stored correctly at (+2 to +8) °C, it is stable till the expiration date stated on the cover. After use the vials must be screwed tightly and stored at (+2 to +8) °C. Avoid contamination, pipette and decantation of standard back to the vial must be excluded.

Precautions in use

The set is specified only for use in vitro. To be handled by entitled and professionally educated person. Serum was tested for presence of HBsAg and antibodies against HIV and HCV with negative results. In spite of the fact, the presence of some pathogenic materials can not be excluded. For this reason the content of the set has to be treated as potentially infectious material.

Reagents of the kit are not classified like dangerous, but contain in low concentration (<0.1%) sodium azide, very toxic substance, dangerous for the environment.

First aid

If swallowed, a half liter of water should be drunk. In case of contact with skin or eyes, rinse them out with drinking water.

The qualified medical aid has to be provided.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging paper waste and containers to recycling.

References

Johnson, A.M.: A new international preparation for proteins in human serum, Arch.Pathol.Lab.Med.117 (1993) 29-31

Date of last revision: 8.12.2014

Produced by Erba Lachema s.r.o.
Karásek 1d, 621 00 Brno, CZ
e-mail: diagnostics@erbalachema.com, www.erbalachema.com**IMULATEST****Multišstandard SET 3****(MSS 3x 1)**

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