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Cat. No. 10010234

Store (+2 to +8) °C

The kit is intended for direct quantitative determination of glycated hemoglobin (HbA1c) in human blood used kit HbA1c Direct 1x0,5, Cat. No. 10010233.

Lot No.:	xxxxxx
Exp. Date:	mm/rrrr

Assigned values according to IFCC

Instrument	Level 1 mmol/mol	Level 2 mmol/mol	Level 3 mmol/mol	Level 4 mmol/mol
Hitachi 717	00,00	00,00	00,00	00,00
Hitachi 917	00,00	00,00	00,00	00,00
Cobas Mira	00,00	00,00	00,00	00,00
on board lyse/Olympus	00,00	00,00	00,00	00,00

Calculation of units: mmol/mol = 10 x % (IFCC)

The following equation can be used for calculation of units IFCC/NGSP:

$$NGSP = (0,915 \times IFCC) + 2,15$$

Characteristics

HbA1c calibrator is lyophilised hemolysate.

It is monospecific calibrator, values are LOT specific.

The setpoint of the calibrators were obtained by assaying representative samples of the entire lots against materials referenced to NGSP. Values for IFCC were recalculated.

Reagents

HbA1c Calibrator 4x 0,5 ml

Composition

Hemolysate was prepared from packed human erythrocytes. Stabilisers for maintain hemoglobin in reduced state for accurate measurement of the HbA1c.

Preparation

1. Open the vial carefully, avoiding any loss of the lyophilised material.
2. Add exactly 0,5 ml of deionised water (inaccurate reconstitution of the control and error in assay technique can caused erroneous results)
3. Close the vial carefully and gently mix for 10 minutes, or until all material has dissolved, avoiding the formation of foam. **Do not shake!**
4. The reconstituted calibrator should be assays in the same manner as the blood specimens including the hemolysate procedure (please, refer to the instruction insert for HbA1c Direct 1x0,5, Cat. No. 10010233):
 - Dispense 1 ml of hemolysis reagent into a tube
 - Add 20 µl of well mixed whole blood and mix
 - Allow to stand for 5 minutes

Stability

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

Stability after reconstitution:

30 days at (+2 to +8) °C in the dark, 3 months at (-20 °C) in the dark
Close vials carefully immediately after use. Protect from the light, heat and contamination. Do not use the product if there is visible evidence of microbial growth of the vial. Improper handling and/or store of the calibrators can affect results.

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Although this product has been tested and found non-reactive for HbsAg and HIV, no known test can offer assurance that products derived from human blood will not transmit diseases. Therefore all human serum products and patients specimens should be handled in the same manner as an infectious agent.

Do not pipette by mouth. Avoid contact with skin and mucous membranes.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0,5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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