**IMULATEST****HbA1c Kontrola****(HbA1c CON 4x0,5)**

Kat. č. 10010235

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Cat. No. 10010235

Store at (+2 to +8) °C

The kit is intended for direct quantitative determination of glycated hemoglobin (HbA1c) in human blood used kit HbA1c Direct 1x0.5, Cat. No. 10010233.

Lot No.:	xxxxxx
Exp. Date:	mm/yyyy

Assigned values according to IFCC

	Mean (mmol/mol)	Range (mmol/mol)
Level 1	00,00	00,00 – 00,00
Level 2	00,00	00,00 – 00,00

Calculation of units: mmol/mol = 10 x %(IFCC)  
 The following equation can be used for calculation of units IFCC/NGSP:  
 NGSP = (0.915 x IFCC) + 2.15

**Characteristics**

HbA1c control is lyophilised hemolysate. It is a monospecific calibrator. It provides 2 levels of HbA1c, one level in the normal range and the other level in the elevated range.

**Reagents**

R1 HbA1c Normal control 2x0.5 ml  
 R2 HbA1c Elevated control 2x0.5 ml

**Composition**

Hemolysate was prepared from packed human erythrocytes. Stabilisers for maintain hemoglobin in reduced state for accurate measurement of the HbA1c.

**Preparation**

1. Open the vial carefully, avoiding any loss of the lyophilised material.
2. Add exactly 0.5 ml of deionised water (inaccurate reconstitution of the control and error in assay technique can caused erroneous results)
3. Close the vial carefully and gently mix for 10 minutes, or until all material has dissolved, avoiding the formation of foam. Do not shake!
4. The reconstituted controls should be assays in the same manner as the blood specimens including the hemolysate procedure (please, refer to the instruction insert for HbA1c Direct 1x0.5, Cat. No. 10010233):
  - Dispense 1 ml of hemolysis reagent into a tube
  - Add 20 µl of well mixed whole blood and mix
  - Allow to stand for 5 minutes

**Stability**

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package. Stability after reconstitution: 30 days at (+2 to +8) °C, in the dark 3 months at (-20) °C, in the dark

Freeze and that only once!

Close vials carefully immediately after use. Protect from the light, heat and contamination.

Do not use the product if there is visible evidence of microbial growth of the vial. Improper handling and/or store of the controls can affect results.

**Health protection**

For in vitro diagnostic use. To be handled by entitled and professionally educated person. Although this product has been tested and found non-reactive for HbsAg and HIV, no known test can offer assurance that products derived from human blood will not transmit diseases. Therefore all human serum products and patients specimens should be handled in the same manner as an infectious agent. Do not pipette by mouth. Avoid contact with skin and mucous membranes.

**First aid**

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

**Waste disposal**

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

**Reference**

13. Trivelli, LA., Ranney, H.M., and Lai, H.T., New Eng. J.Med. 284, 353 (1971)
14. Gonen, B., and Rubenstein, AH., Diabetologia 15, 1 (1978)
15. Gabbay, K.H., Hasty, K., Breslow, J.L., Ellison, R.C., Bunn, H.F., and Gallop, P.M., J. Clin. Endocrinol. Metab. 44, 859 (1977)
16. Bates, H.M., Lab. Mang., Vol 16 (Jan 1978)

Date of last revision: 12. 8. 2014

Выпускает Erba Lachema s. r. o., Karásek 1d, 621 00 Brno, CZ  
 e-mail: diagnostics@erbalachema.com, www.erbalachema.com

Produce Erba Lachema s. r. o., Karásek 1d, 621 00 Brno, CZ  
 e-mail: diagnostics@erbalachema.com, www.erbalachema.com

IM/PI/51/14/S