IMULATEST**CYSTATIN C Kalibrátor (CYS C CAL 1x2)**

Kat. č. 10010249

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Cat. No. 10010249

Store at (+2 to +8) °C

The kit is intended for immunoturbidimetric quantitative determination of Cystatin C in human serum.

Batch	xxxxxx
Expiration	mm/yyyy

Certified values (mg/l)

Parameter	Level
Cystatin C	0.00

Characteristics

CYSTATIN C calibrator is a synthetic liquid calibration material, that contains human Cystatin C. It is used for calibration during determination of Cystatin C in human serum. It contains sodium azide (0,095%).

Way of use and stability

Cystatin C calibrator is ready for use. It is treated as calibrator according to the instruction included to the diagnostic set. Prior use, allow Cystatin C calibrators to equilibrate to room temperature (+22 to +28)°C.

If the kit is stored correctly at (+2 to +8)°C, it is stable till the expiration date stated on the cover. After use the vials must be screwed tightly and stored at (+2 to +8) °C.

After the first opening, the stability is 4 weeks at (+2 to +8) °C.

Do not freeze!

Cystatin C calibrators should not be allowed for repeated long periods of time (up to 6 hours) at room temperature.

Precautions in use

The set is specified only for use in vitro. To be handled by entitled and professionally educated person. Serum was tested for presence of HBsAg and antibodies against HIV and HCV with negative results. In spite of the fact, the presence of some pathogenic materials can not be excluded. For this reason the content of the set has to be treated as potentially infectious material.

Reagents of the kit are not classified like dangerous, but contain in low concentration (0,095%) sodium azide, very toxic substance, dangerous for the environment.

First aid

If swallowed, a half liter of water should be drunk. In case of contact with skin or eyes, rinse them out with drinking water.

The qualified medical aid has to be provided.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and containers to recycling.

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Tendstad O, Roald AB, Grubb A, Auckland K, Scan.J.Clin.Lab.Invest. 1996, 56:409-14