IMULATEST**ASLO Kalibrátor****(ASLO CAL 1x 1)**

Kat. č. 10010252

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Literatúra

Spaun, J. et al., International standard for anti-streptolysin-O. Bull. WHO 24 (1961) 271-279.

Datum poslednej revízie: 12. 8. 2014

Vyrába Erba Lachema s.r.o.
Karásek 1d, 621 00 Brno, CZ
e-mail: diagnostics@erbalachema.com, www.erbalachema.com

IMULATEST**ASLO Calibrator (ASLO CAL 1x 1)**

Cat. No. 10010252 Store at (+2 to +8) °C

The kit is intended for immunoturbidimetric quantitative determination of antistreptolysin O in human serum .

Batch	xxxxxx
Expiration	mm/yyyy

Certified values

Parameter	Concentration	Uncertainty	Unit
ASLO	000	00	IU/ml

Referential material, where the values are related to:
International standard SSI
(1st International Standard for Antistreptolysin O)

Characteristics

ASLO calibrator is a synthetic liquid calibration material, that contains human Antistreptolysin O. It is used for calibration during determination of Antistreptolysin O in human serum. It contains sodium azide (0,075%).

Way of use and stability

The preparation is ready for use. It is treated as calibrator according to the instruction included to the diagnostic set.

If the kit is stored correctly at (+2 to +8) °C, it is stable till the expiration date stated on the cover. After use the vials must be screwed tightly and stored at (+2 to +8) °C. Avoid contamination, pipette and decantation of calibrator back to the vial must be excluded.

Precautions in use

The set is specified only for use in vitro. To be handled by entitled and professionally educated person. Serum was tested for presence of HBsAg and antibodies against HIV and HCV with negative results. In spite of the fact, the presence of some pathogenic materials can not be excluded. For this reason the content of the set has to be treated as potentially infectious material.

First aid

If swallowed, a half liter of water should be drunk. In case of contact with skin or eyes, rinse them out with drinking water. The qualified medical aid has to be provided.

Waste disposal

All tested samples should be treated as potentially infectious and with the contingent rest of the reagents should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging paper waste and containers to recycling.

References

Spaun, J. et al., International standard for anti-streptolysin-O. Bull. WHO 24 (1961) 271-279.

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