LYONORM LIPID HUM N

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(EN)

Characteristics

LYONORM LIPID HUM N is a lyophilized control serum based on human matrix with analyte concentrations in the normal up to slightly increased range.

Use

The control serum LYONORM LIPID HUM N is intended especially for control of accuracy and precision of the determination of HDL and LDL cholesterol and other analytes (cholesterol, triglycerides, apolipoproteins) by HDL C DIRECT 240 (Cat. No. BLT00030) and LDL C DIRECT 80 (Cat. No. BLT00043) and 240 (Cat. No. BLT00042) tests.

Values and ranges

The actual values and ranges of concentrations of individual analytes are specified in the certificate of analysis attached. The lot number of the bottle must match with the number given in the certificate.

Reconstitution of lyophilisate

Open the bottle very carefully avoiding any loss of lyophilized material. Add exactly 5.00 ml of redistilled water. Close the bottle tightly, swirl gently and allow to stand for 30 min in the dark. Let the content to dissolve completely avoiding any formation of the foam. Do not shake the bottle.

Stability and storage

When stored at a dark place, the LYONORM LIPID HUM N is stable in the lyophilized form:

at temperature of (+2 to +8) °C till the date of expiration stated on the bottle.

The reconstituted serum is stable for at least:

- 1 day at (+15 to +20) °C in dark
- 7 days at (+2 to +8) °C in dark
- 1 month at (-20) °C in dark, when frozen once

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

The reagents of the set contain biological material, which is potentially infectious.

The control sera have been prepared from normal human serum, which was tested for HIV and HBsAg antibodies by ELISA or RIA methods with negative results.

The negative results of tests for HIV antibodies in the case of lyophilized products were confirmed by immunofluorescent method performed by the authorized institution.

However, since the risk of infection can never be excluded, the product is to be handled like potentially infectious human material, similar to sera of patients.

During the work all principles of general safety are to be observed. Above all observe the principles of personal hygiene, use personal protection aids (clothes, shoes, etc.), do not eat, drink and smoke. After finishing the work and before a meal wash your hands with soap and warm water and apply a suitable regenerating cream.

When the bottle is broken, the contaminated place must be disinfected with Persteril or Chloramin.

Note

LYONORM LIPID HUM N is control serum for in vitro diagnostic use. Not to be used internally in humans or animals. Do not use this product past the expiration date stated on the bottle. Do not use this product for any purpose other than described herein.

Waste disposal

The kit should be held as a potentially infectious material. The eventual rest of the reagent must be disposed according to inner regulations for dangerous waste, in compliance with local and national regulations. Paper packing are handed over for recycling or discarded as communal waste.

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