LIPOPROTEIN (a)

(IT-LP (a) 1 x 2)

Kat. č. 10003242

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Cat. No. 10003242

Store at (+2 to +8) °C

Reagent kit for immunoturbidimetric determination of lipoprotein (a) in human serum and plasma.

Principle

The human serum (plasma) sample reacts upon a specific antibody for human Lp (a) and the turbidity induced by the formation of immune complexes is recorded at 340 nm. The turbidity measured is directly proportional to the Lp (a) sample concentration.

Reagents

R1 Buffer

(1 x 30 ml)

Imidazole buffer pH 7.0	35 mmol/l
PEG	40 g/l
Sodium chloride	0.15 mol/l
Sodium azide	1 g/l

R2 Antiserum

(1 x 2 ml)

Sterile delipidated goat serum against human Lp (a) titre +/- 4 mg/ml

Reaction mixture

Imidazole buffer pH 7.0	31.9 mmol/l
PEG	36.5 g/l
Sodium chloride	136.8 mmol/l
Sodium azide	< 1 g/l

Preparation and stability of working reagents

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

Reagents R1 and R2 are liquid ready to use. After opening, reagent R1 is stable until expiry date at (+15 to +25) °C and reagent R2 is stable until expiry date at (+2 to +8) °C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

Samples

Serum or plasma (do not use heparin or EDTA plasma), samples should be fresh and not hemolysed.

Prior the analysis, the serum can be stored a week at (+2 to +8) °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

Calibration

For calibration use Lp (a) standard (Cat. No. 10003243).

For photometers working with a 6 points calibration curves, prepare a range of 6 reference samples by successive 1:1 dilutions of the Lp (a) standard by saline. If, for example, the concentration of Lp (a) in the standard is 2500 mg/l, prepare the following concentrations: 75,56-151,25-312,5-625-1250-2500 mg/l (dilution in ratio 1:32,1:16,1:8,1:4, 1: 2, undiluted standard). Plot the calibration curve.

Dilute samples with Lp (a) concentration higher than the highest standard by saline 10x and repeat analysis. The results multiply by dilution.

Quality control

For quality control use Bio-Rad Liquicheck™ Immunology Control or other quality control sera with Lp (a) concentration values determined by immunoturbidimetric method.

Procedure

Wavelength	340 nm
Cuvette	1 cm
Temperature	37 °C
Serum/reaction mixture ratio	1/68.5
Reagents and sample volume can be modified, by respecting reagents/sample volume ratio	

	Sample	Standard	Reagent blank
Standard	–	0.003 ml	–
Sample	0.003 ml	–	–
Distilled water	–	–	0.003 ml
Reagent 1	0.250 ml	0.250 ml	0.250 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A1.			
Reagent 2	0.025 ml	0.025 ml	0.025 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A2.			

Applications for automatic analysers can be supplied for request.

Calculation

Calculate the difference A2-A1 for standards and samples.

$$LP(a) \text{ (mg/l)} = \frac{(A_2 - A_1)_{\text{sample}} - (A_2 - A_1)_{\text{reagent blank}}}{(A_2 - A_1)_{\text{st}} - (A_2 - A_1)_{\text{reagent blank}}} \times C_{\text{st}}$$

C_{st} ...concentration of standard

Reference values

fS Lp (a) (mg/l) < 300

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

Performance characteristics

Limit of detection: 30 mg/l

Hook effect: from 2 800 mg/l

Specificity: the antiserum shows the unique line against the whole human serum in the immunoelectrophoresis and in the passive immunodiffusion.

Repeatability: < 6 %, **Reproducibility:** < 9 %

Interferences: Triglycerides to < 11.3 mmol/l, hemoglobin to < 200 µmol/l, bilirubin to < 500 µmol/l do not interfere with the measurements (< 10% or 2 SD).

Calibration stability: on the most often automatic analysers it is < 3 weeks.

Health protection

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.

Reagents of the kit are not classified as dangerous however, they contain in low concentration a toxic sodium azide with exposure limit. Reagent 1 also contains 0.2% imidazole which may damage an unborn child.

EUH210 Safety data sheet available on request.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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