C-REAKTIVNÍ PROTEIN READY TO USE

(IT-CRP DIL 1x15)

Kat.č. 10004971

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Literatúra

- Turner, M. W. & Hulme, B. The Plasma proteins: An introduction. Pittman Medical & Scientific, London, 1971.
- The Acute phase response, Abstracts of the Biochemical Society meeting, Belfast Sept. 14-17, The Biochemical Society, London 1993
- Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Borquez de Larrea, L. et. al., Determination of C-reaktívni protein by an improved turbidimetric assay on Boehringer Mannheim/Hitachi analysis system, Klin.Lab. 39 (1993) 55-62
- Zima, T.: Laboratórna diagnostika, Galén (2002)

Dátum poslednej revízie: 25. 10. 2013

C-REACTIVE PROTEIN READY TO USE

(IT-CRP DIL 1x15)

Cat. No. 10004971

Store at (+2 to +8) °C

Reagent kit for immunoturbidimetric determination of C-reactive protein in human serum and plasma.

Principle

The human serum (plasma) sample reacts upon a specific antibody for human CRP and the turbidity induced by the formation of immune complexes is recorded at 340 nm. The turbidity measured is directly proportional to the CRP sample concentration.

Reagents

R1 Buffer	(1 x 125 ml)
Phosphate buffer pH 7.4	35 mmol/l
PEG	40 g/l
Sodium chloride	150 mmol/l
Sodium azide	0.9 g/l

R2 Antiserum	(1x15 ml)
Sterile delipidated goat serum against human CRP	titre +/- 1 g/l
HEPES buffer pH 7.4	50 mmol/l
EDTA	9 mmol/l
Sodium azide	0.9 g/l

Reaction mixture

Phosphate buffer pH 7.4	28.2 mmol/l
HEPES buffer pH 7.4	6.5 mmol/l
PEG	32.2 g/l
Sodium chloride	121 mmol/l
EDTA	1.16 mmol/l
Sodium azide	0.7 g/l

Preparation and stability of working reagents

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

Reagents R1 and R2 are liquid ready to use. After opening, reagent R1 is stable until expiry date at (+15 to +25) °C and reagent R2 is stable until expiry date at (+2 to +8) °C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

Samples

Serum or plasma, samples should be fresh and not hemolysed. Prior the analysis, the serum can be stored a week at (+2 to +8) °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

Calibration

For calibration use CRP standard (Cat. No. 10003237). For photometers working with a 6 points calibration curves, prepare a range of 6 reference samples by successive 1:1 dilutions of the CRP standard by saline. If, for example, the concentration of CRP in the standard is 200 mg/l, prepare the following concentrations: 6.25-12.5-25-50-100-200 mg/l (dilution in ratio 1:32,1:16,1:8, 1:4, 1:2, undiluted standard). Plot the calibration curve. Dilute samples with CRP concentration higher than the highest standard by saline 10x and repeat analysis. The result multiply by dilution.

Quality control

For quality control use Bio-Rad Liquicheck™ Immunology Control or other quality control sera with CRP concentration values determined by immunoturbidimetric method.

Procedure

Wavelength	340 nm
Cuvette	1 cm
Temperature	37 °C
Serum/reaction mixture ratio	1/15.5

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

	Sample	Standard	Reagent blank
Standard	–	0.020 ml	–
Sample	0.020 ml	–	–
Distilled water	–	–	0.020 ml
Reagent 1	0.250 ml	0.250 ml	0.250 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A1..			
Reagent 2	0.040 ml	0.040 ml	0.040 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A2.			

Applications for automatic analysers can be supplied for request.

Calculation

Calculate the difference A2-A1 for standards and samples.

$$\text{CRP (mg/l)} = \frac{(A_2 - A_1)_{\text{sample}} - (A_2 - A_1)_{\text{reagent blank}}}{(A_2 - A_1)_{\text{st}} - (A_2 - A_1)_{\text{reagent blank}}} \times C_{\text{st}}$$

C_{st} ...concentration of standard

Reference values

fS C-reactive protein (mg/l) < 5

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

Performance characteristics

Limit of detection: 2.4 mg/l

Hook effect: from 950 mg/l

Specificity: the antiserum show the unique line against the whole human serum in the immunoelectrophoresis and in the passive immunodiffusion.

Precision (at 37 °C)

INTRA-ASSAY	Mean	SD	CV
n = 10	(mg/l)	(mg/l)	(%)
Sample 1	12.7	0.67	5.3
Sample 2	22.1	0.74	3.3
INTER-ASSAY	Mean	SD	CV
n = 10	(mg/l)	(mg/l)	(%)
Sample 1	13.4	0.49	3.7
Sample 2	22.8	0.91	4.0

Comparison: with commercial available test (immunoturbidimetry). Results: n = 54, r = 0.996, y = 1.02x + 1.8 mg/l

Interferences: Triglycerides to < 11.3 mmol/l, hemoglobin to < 200 mmol/l, bilirubin to < 100 mmol/l do not interfere with the measurements (< 10% or 2 SD).

Calibration stability: on the most often automatic analysers it is > 4 weeks.

Health protection

For in vitro diagnostic use. To be handled by entitled and professionals educated person.

Reagent of the kit are not classified like dangerous.

First aid

In case of an accidental ingestion wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

References

- Turner, M. W. & Hulme, B. The Plasma proteins: An introduction. Pittman Medical & Scientific, London, 1971.
- The Acute phase response, Abstracts of the Biochemical Society meeting, Belfast Sept. 14-17, The Biochemical Society, London 1993
- Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Borque de Larrea, L. et. al., Determination of C-reaktivní protein by an improved turbidimetric assay on Boehringer Mannheim/Hitachi analysis system, Klin. Lab. 39 (1993) 55-62
- Zima, T.: Laboratorní diagnostika, Galén (2002)

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1. Turner, M. W. & Hulme, B. The plasma proteins: An introduction. Pittman Medical & Scientific, London, 1971.
2. The Acute phase response. Abstracts of the Biochemical Society meeting, Belfast Sept. 14-17. The Biochemical Society, London, 1993.
3. Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
4. Borque de Larrea, L. et al. Determination of C-reactive protein by an improved turbidimetric assay on Boehringer Mannheim/Hitachi analysis systems. Klin. Lab. 39 (1993) 55-62.

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