IMUNOGLOBULIN M READY TO USE

(IT-IgM DIL 1x15)

Kat. č. 10006878

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Cat. No. 10006878

Store at (+2 to +8) °C

Reagent kit for immunoturbidimetric determination of immunoglobulin M in human serum and plasma.

Principle

The human serum (plasma) sample reacts upon a specific antibody for human IgM and the turbidity induced by the formation of immune complexes is recorded at 340 nm. The turbidity measured is directly proportional to the IgM sample concentration.

Reagents

R1 Buffer

Tricine buffer pH 8	(1x125 ml)	35 mmol/l
PEG		40 g/l
Sodium chloride		154 mmol/l
Sodium azide		0.9 g/l

R2 Antiserum

Sterile delipidated goat serum against human IgM	(1x15 ml)	titre +/- 6 g/l
HEPES buffer pH 7.4		50 mmol/l
EDTA		9 mmol/l
Sodium azide		0.9 g/l

Reaction mixture

Tricine buffer pH 8	28.88 mmol/l
HEPES buffer pH 7.4	8.25 mmol/l
PEG	33 g/l
Sodium chloride	127.1 mmol/l
EDTA	1.49 mmol/l
Sodium azide	0.89 g/l

Preparation and stability of working reagents

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

Reagents R1 and R2 are liquid ready to use. After opening, reagent R1 is stable until expiry date at (+15 to +25) °C and reagent R2 is stable until expiry date at (+2 to +8) °C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

Samples

Serum or plasma (heparin, NaF, EDTA, sodium citrate), samples should be fresh and not hemolysed.

Prior the analysis, the serum can be stored a week at (+2 to +8) °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

Calibration

For calibration use Multistandard high (Cat. No. 10003233) or Multistandard set 3 (Cat. No. 10003238).

For photometers working with a 6 points calibration curves, prepare a range of 6 reference samples by successive 1:1 dilutions of the Multistandard high by saline. If, for example, the concentration of IgM Multistandard high is 3.0 g/l, prepare the following concentrations: 0.10-0.19-0.38-0.75-1.5-3.0 g/l (dilution in ratio 1:32, 1:16, 1:8, 1:4, 1:2, undiluted standard).

Plot the calibration curve.

Dilute samples with IgM concentration higher than the highest standard by saline 10x and repeat analysis. The result multiply by dilution.

Quality control

For quality control use Bio-Rad Liquicheck™ Immunology Control or other quality control sera with IgM concentration values determined by immunoturbidimetric method.

Procedure

Wavelength	340 nm
Cuvette	1 cm
Temperature	37 °C
Serum/reaction mixture ratio	1/101

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

	Sample	Standard	Reagent blank
Standard	–	0.003 ml	–
Sample	0.003 ml	–	–
Distilled water	–	–	0.003 ml
Reagent 1	0.250 ml	0.250 ml	0.250 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A1.			
Reagent 2	0.050 ml	0.050 ml	0.050 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A2.			

Applications for automatic analysers can be supplied for request.

Calculation

Calculate the difference A2-A1 for standards and samples.

$$\text{IgM(g/l)} = \frac{(A_2 - A_1)_{\text{sample}} - (A_2 - A_1)_{\text{reagent blank}}}{(A_2 - A_1)_{\text{standard}} - (A_2 - A_1)_{\text{reagent blank}}} \times C_{\text{st}}$$

C_{st} ... concentration of standard

Reference values

fS IgM (g/l)

adults	0.4 – 2.1
1 year to 12 years children	0.4 – 1.5
1 month to 12 months children	0.1 – 1.0

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

Performance characteristics

Limit of detection: 0.005 g/l

Hook effect: from 20.0 g/l

Specificity: the antiserum show the unique line against the whole human serum in the immunoelectrophoresis and in the passive immunodiffusion.

Precision (at 37 °C)

INTRA-ASSAY	Mean (g/l)	SD (g/l)	CV (%)
n = 10			
Sample 1	0.62	0.009	1.5
Sample 2	1.24	0.011	0.9
INTER-ASSAY	Mean (g/l)	SD (g/l)	CV (%)
n = 10			
Sample 1	0.57	0.028	5.0
Sample 2	1.24	0.067	5.4

Comparison: with commercial available test (nephelometry).

Results: n = 31, r = 0.980, y = 1.1x + 0.018 g/l

Interferences: Triglycerides to < 11.3 mmol/l, hemoglobin to < 200 µmol/l, bilirubin to < 500 µmol/l do not interfere with the measurements (< 10% or 2 SD).

Calibration stability: on the most often automatic analysers it is > 4 weeks.

Health protection

For in vitro diagnostic use. To be handled by entitled and professionals educated person.

Reagent of the kit are not classified like dangerous.

Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.

Products containing sodium azide have to be handled with care, avoid ingestion and contact with skin and mucous membrane.

First aid

In case of an accidental ingestion wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

References

- Bergström, K. & Lefvert, A. K.: Scand. J. Clin. Lab. Invest. 40 (1980) 637
- Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
- Roitt, I.: Essential Immunology, Blackwell, Oxford (1991)
- Zima, T.: Laboratorní diagnostika, Galén (2002)

USED SYMBOLS

Catalogue Number

Manufacturer

See Instruction for Use

Lot Number

Storage Temperature

Expiry date

In vitro Diagnostics

Content

CE Mark - Device comply with the Directive 98/79/EC

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