## TRANSFERIN READY TO USE

(IT-TRF DIL 1x15)

Kat. č. 10007748

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Cat. No. 10007748

Store at (+2 to +8) °C

Reagent kit for immunoturbidimetric determination of transferrin in human serum and plasma.

### Principle

The human serum (plasma) sample reacts upon a specific antibody for human transferrin and the turbidity induced by the formation of immune complexes is recorded at 340 nm. The turbidity measured is directly proportional to the transferrin sample concentration.

### Reagents

<b>R1 Buffer</b>	<b>(1x125 ml)</b>
Imidazol buffer pH 7	35 mmol/l
PEG	40 g/l
Sodium chloride	151 mmol/l
Sodium azide	0.9 g/l
<b>R2 Antisérúm</b>	<b>(1x15 ml)</b>
Sterile delipidated goat serum against human transferrin	titr +/- 2 g/l
HEPES buffer pH 7.4	50 mmol/l
EDTA	9 mmol/l
Sodium azide	0.9 g/l

### Reaction mixture

Imidazol buffer pH 7	27 mmol/l
HEPES buffer pH 7.4	11 mmol/l
PEG	31 g/l
Sodium chloride	118 mmol/l
EDTA	2 mmol/l
Sodium azide	0.9 g/l

### Preparation and stability of working reagents

If stored at (+2 to +8) °C, reagents are stable until expiry date, that is stated on the package.

Reagents R1 and R2 are liquid ready to use. After opening, reagent R1 is stable until expiry date at (+15 to +25) °C and reagent R2 is stable until expiry date at (+2 to +8) °C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

### Samples

Serum or plasma (heparin, NaF, EDTA, sodium citrate), samples should be fresh and not hemolysed.

Prior the analysis, the serum can be stored a week at (+2 to +8) °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

### Calibration

For calibration use Multistandard high (Cat. No. 10003233) or Multistandard set 3 (Cat. No. 10003238).

For photometers working with a 5 points calibration curves, prepare a range of 5 reference samples by successive 1:1 dilutions of the Multistandard high by saline. If, for example, the concentration of transferrin in Multistandard high is 7.0 g/l, prepare the following concentrations: 0.44-0.88-1.75-3.5-7.0 (dilution in ratio 1:16, 1:8, 1:4, 1:2, undiluted standard).

Plot the calibration curve.

Dilute samples with transferrin concentration higher than the highest standard by saline 10x and repeat analysis. The results multiply by dilution.

### Quality control

For quality control use Bio-Rad Liquicheck™ Immunology Control or other quality control sera with transferrin concentration values determined by immunoturbidimetric method.

### Procedure

Wavelength	340 nm
Cuvette	1 cm
Temperature	37 °C
Serum/reaction mixture ratio	1/161

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

### SYMBOLS USED ON LABEL

REF	Catalogue Number	Manufacturer	See Instruction for Use	LOT	Lot Number
Storage Temperature	Expiry date	IVD	In vitro Diagnostics	CONT	Content

CE Mark - Device comply with the Directive 98/79/EC

	Sample	Standard	Reagent blank
Standard	—	0.002 ml	—
Sample	0.002 ml	—	—
Distilled water	—	—	0.002 ml
Reagent 1	0.250 ml	0.250 ml	0.250 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A1.			
Reagent 2	0.070 ml	0.070 ml	0.070 ml
Mix and incubate for 5 min. at temperature of 37 °C. Then measure the absorbance A2.			

Applications for automatic analysers can be supplied for request.

### Calculation

Calculate the difference A2-A1 for standards and samples.

$$\text{TRF (g/l)} = \frac{(A_2 - A_1)_{\text{sample}} - (A_2 - A_1)_{\text{reagent blank}}}{(A_2 - A_1)_{\text{st}} - (A_2 - A_1)_{\text{reagent blank}}} \times C_{\text{st}}$$

C<sub>ST</sub> ... concentration of standard

### Reference values

fS transferrin (g/l) 1,69 – 3,09

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

### Performance characteristics

**Limit of detection:** 0,091 g/l

**Hook effect:** from 16,8 g/l

**Specificity:** the antiserum shows the unique line against the whole human serum in the immunoelectrophoresis and in the passive immunodiffusion.

### Precision (at 37 °C)

INTRA-ASSAY	Mean	SD	CV
n = 10	(g/l)	(g/l)	(%)
Sample 1	2.44	0.04	1.6
Sample 2	3.74	0.06	1.5
INTER-ASSAY	Mean	SD	CV
n = 10	(g/l)	(g/l)	(%)
Sample 1	2.42	0.08	3.4
Sample 2	3.89	0.08	2.1

**Comparison:** with commercial available test (nephelometry).

**Results:** n = 36, r = 0,941, y = 0,96 x + 0,3 g/l

**Interferences:** Triglycerides to < 22,6 mmol/l, hemoglobin to < 160 µmol/l, bilirubin to < 500 µmol/l do not interfere with the measurements (< 10% or 2 SD).

Calibration stability: on the most often automatic analysers it is >4 weeks.

### Health protection

For in vitro diagnostic use. To be handled by entitled and professionals educated person.

Reagents of the kit are not classified like dangerous.

Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.

Products containing sodium azide have to be handled with care, avoid ingestion and contact with skin and mucous membrane.

### First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

### Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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