IMULATEST

CZ

ALBUMIN (V MOČI) READY TO USE

(IT-ALB DIL 1x15)

Kat. č. 10007760

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Dátum poslednej revízie: 4.12. 2015

IVD

ALBUMIN (IN URINE) READY TO USE (IT-ALB DIL 1x15)

Cat. No. 10007760

Store at (+2 to +8)°C

Reagent kit for immunoturbidimetric determination of albumin in urine (microalbuminuria).

Principle

The human serum (plasma) sample reacts upon a specific antibody for human albumin and the turbidity induced by the formation of immune complexes is recorded at 340 nm. The turbidity measured is directly proportional to the albumin sample concentration.

Reagents

R1 Buffer	(1x125 ml)
Imidazol buffer pH 7	35 mmol/l
PEG	40 g/l
Sodium chloride	151 mmol/l
Sodium azide	0.9 g/l
R2 Antisérum	(1x15 ml)
Sterile delipidated goat serum against human albumin	titr +/- 1.5 g/l
HEPES buffer pH 7.4	50 mmol/l
EDTA	9 mmol/l
Sodium azide	0.9 g/l
Reaction mixture	
Imidazol buffer pH 7	29 mmol/l
HEPES buffer pH 7.4	8 mmol/l
PEG	33 g/l
Sodium chloride	126 mmol/l
EDTA	1.5 mmol/l
Sodium azide	0.9 g/l

Preparation and stability of working reagents

If stored at (+2 to +8)°C, reagents are stable until expiry date, that is stated on the package.

Reagents R1 and R2 are liquid ready to use. After opening, reagent R1 is stable until expiry date at (+15 to +25)°C and reagent R2 is stable until expiry date at (+2 to +8)°C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

Samples

Urine, samples should be fresh, without preservative. Prior the analysis, the sample can be stored a week at (+2 to +8)°C or 3 months at -20°C. Avoid repeated freezing and thawing.

Calibration

For calibration use MAU standard (Cat. No. 10003230). For photometers working with a 6 points calibration curves, prepare a range of 6 reference samples by successive 1:1 dilutions of the MAU standard by saline. If, for example, the concentration of transferrin in MAU standard is 250 mg/l, prepare the following concentrations: 7.8-15.6-31.3-62.5-125-250 mg/l (dilution in ratio 1:32, 1:16, 1:8, 1:4, 1:2, undiluted standard). Plot the calibration curve. Dilute samples with albumin concentration higher than the highest standard by saline 10x and repeat analysis. The results multiply by dilution.

Quality control

For quality control use Bio-Rad Liquicheck™ Immunology Control or other quality control sera with albumin concentration values determined by immunoturbidimetric method.

Procedure

Wavelength	340 nm
Cuvette	1 cm
Temperature	37°C
Sample/reaction mixture ratio	1/31

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio

	-	0.010 ml	-
Standard	0.010 ml	-	-
Sample	-	-	0.010 ml
Distilled water	-	0.010 ml	-
Reagent 1	0.250 ml	0.250 ml	0.250 ml
Mix and incubate for 5 min. at temperature of 37°C. Then measure the absorbance A1.			
Reagent 2	0.050 ml	0.050 ml	0.050 ml
Mix and incubate for 5 min. at temperature of 37°C. Then measure the absorbance A2.			

Applications for automatic analysers can be supplied for request.

Calculation

Calculate the difference A2-A1 for standards and samples.

$$\text{ALBUMIN (mg/l)} = \frac{(A_2 - A_1)_{st} - (A_2 - A_1)_{\text{reagentblank}}}{(A_2 - A_1)_{st} - (A_2 - A_1)_{\text{reagentblank}}} \times C_{st}$$

C_{ST} ... concentration of standard

Reference values

fU albumin (mg/l) < 20
The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

Performance characteristics

Limit of detection: 3.6 mg/l

Hook effect: from 1200 mg/l

Specificity: the antiserum shows the unique line against the whole human serum in the immunoelectrophoresis and in the passive immunodiffusion.

Precision (at 37°C)

INTRA-ASSAY	Mean (mg/l)	SD (mg/l)	CV (%)
n = 10			
Sample 1	6.57	0.48	7.3
Sample 2	31.08	0.90	2.9
INTER-ASSAY	Mean (mg/l)	SD (mg/l)	CV (%)
n = 10			
Sample 1	6.66	0.34	5.1
Sample 2	31.75	1.15	3.6

Comparison: with commercial available test (nephelometry).

Results: n = 25, r = 0.994, y = 0.95 x - 1.97 mg/l

Interferences: IgG to <0.03 g/l, bilirubin to < 100 µmol/l do not interfere with the measurements (< 10% or 2 SD).

Calibration stability: on the most often automatic analysers it is >4 weeks.

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous, but contain in low concentration (<0.1%) sodium azide, very toxic substance, dangerous for the environment.

Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

Literatura

- Bergström, K. & Lefvert, A. K.: Scand. J. Clin. Lab. Invest. 40 (1980) 637
- Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1994.
- Roitt, I.: Essential Immunology, Blackwell, Oxford (1991)
- Zima, T.: Laboratorní diagnostika, Galén (2002)

USED SYMBOLS

Catalogue Number

Manufacturer

See Instruction for Use

Lot Number

CE Mark - Device comply with the Directive 98/79/EC

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