IVD

HbA1c Hemolyzační činidlo (HbA1c HEM 500)

Kat. č. 10010236

Skladovat (+2 až +25)°C

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Cat. No. 10010236

Store at (+2 to +25)°C

The kit is intended for the preparation of samples for direct quantitative determination of glycated hemoglobin (HbA1c) in human blood used kit HbA1c Direct 1x0.5, Cat. No. 10010233.

Reagents

R1 Hemolysis reagent	5 x 100 ml
Sodium azide	0.5 g/l

Preparation and stability of working reagents

Reagent is liquid, ready to use.

If stored at (+2 to +25)°C in the dark, reagents are stable until expiry date, that is stated on the package. Close immediately after use.

Please, follow the procedure described at the instruction insert for HbA1c Direct 1x0.5, Cat. No. 10010233.

Hemolysate preparation

1. Dispense 1 ml of hemolysis reagent into a tube
2. Add 20 µl of well mixed whole blood and mix
3. Allow to stand for 5 minutes
4. Follow the same procedure with sets of calibrators and controls

Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous but contain less than 0.1% of sodium azide - classified as very toxic and dangerous substance to the environment.

Products containing sodium azide have to be handled with care, avoid ingestion and contact with skin and mucous membrane.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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CE Mark -
Device comply with
the Directive 98/79/EC

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