IVD

CYSTATIN C

(CYS C 1x5)

Kat. č. 10010248

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Cat. No. 10010248

Store at (+2 to +8) °C

The liquid reagents for immunoturbidimetric quantitative determination of Cystatin C in human serum and plasma.

Principle

The assay of Cystatin C is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immune complexes.

Reagents

R1 Buffer (1x25 ml)
 Tris buffer 1.2 %
 MIT (metylisothiazolinone) 0.03 %

R2 Latex (1x5 ml)
 Suspension of anti-human Cystatin C antibodies coated cystatin C latex particles 0.12 %
 Glycine buffer 1.26 %
 Sodium azide 0.095 %

Preparation and stability of working reagents

If stored at (+2 to +8) °C in the dark, reagents are stable until expiry date, that is stated on the package. Do not freeze the reagents! Reagents R1 and R2 are liquid ready to use.

After first opening, reagent is stable until expiry date at (+2 to +8) °C if stored at appropriate conditions, closed carefully and without any contamination.

Samples

Human serum or plasma (EDTA, heparin)

Prior to the analysis, the serum can be stored 2 days at (+2 to +8) °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

Calibration

It is recommended nonlinear five-point calibration using Cystatin C Calibrator, Cat. No. 10010249, from which the successive dilution in the ratio of 1:1 are prepared in total 5 standards.

Quality control

For quality control, it is recommended to use:

Cystatin C Control, Cat. No. 10010250 (it contains 3x2 levels of Cystatin C, 6x2 ml)

Procedure

Wavelength 546 nm
 Cuvette 1 cm
 Temperature +37 °C
 Sample/reaction mixture ratio 1/94.3
 Reagents and sample volume can be modified provided reagents/sample is respected.

Pipetting schema

Reagent R1 230 µl
 Sample 3 µl
 Pre-incubation 5 min
Reagent R2 50 µl
 Read turbidance A1
 Incubation 5 min
 Read turbidance A2

Hitachi 911 – Automatic analyser setting

TEST CODE	[CYS C]
ASSAY CODE	[2-POINT END]:[10]-[0]
ASSAY POINT	[18]-[31]-[0]]-[0]
SAMPLE VOLUME	[3]
R1 VOLUME	[230]
R3 VOLUME	[50]
WAVELENGTH [2nd/prim]	[800][546]
STD (1) CONC.-POS.	[0.00*][1]
STD (2) CONC.-POS.	[**][2]

STD (3) CONC.-POS.	[**][3]
STD (4) CONC.-POS.	[**][4]
STD (5) CONC. POS.	[**][5]
STD (6) CONC.-POS.	[**][6]
SD LIMIT	[999]
DUPLICATE LIMIT	[1000]
SENSITIVITY LIMIT	[0]
ABS.LIMIT (INC/DEC)	[32000][INCREASE]
PROZONE LIMIT	[-]-[-]
EXPECTED VALUE	[-]-[-]
PANIC VALUE	[-]-[-]
INSTRUMENT FACTOR	[1.00]

* = CALIBRATION 0.0: saline (0,9% NaCl in distilled water)

** = Cystatin C calibrator values

Applications for automatic analysers will be supplied on request.

Reference values

fS (cysC) 0.8 – 1.2 mg/l

The range of reference values is only approximate; it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

Performance data

Limit of detection: 0.09 mg/l

Low limit of quantification: 0.27 mg/l

Hook effect: >30 mg/l

Analytical range: 0.27 – 10.0 mg/l

Precision (at 37 °C)

INTRA-ASSAY n = 10	Mean (mg/l)	SD (mg/l)	CV (%)
Sample 1	1,019	0,028	3,53
Sample 2	2,956	0,050	2,44
INTER-ASSAY n = 10	Mean (g/l)	SD (g/l)	CV (%)
Sample 1	1,019	0,036	2,74
Sample 2	2,956	0,072	1,69

Comparison: with commercial available method. Linear regression:

n = 30, r = 0.995, y = 1.161 x – 0.1519 mg/l

Interference

Following substances do not interfere:

bilirubin up to 410 µmol/l, haemoglobin up to 5 g/l, triglycerides up to 24 mmol/l and ascorbic acid up to 2.8 mmol/l.

Health protection

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified as dangerous. Reagent R2 contains less than 0.1% sodium azide - classified as very toxic and dangerous substance to the environment.

First aid

In case of an accidental ingestion, wash up the mouth and drink about 0,5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

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Tendstad O, Roald AB, Grubb A, Auckland K, Scan.J.Clin.Lab.Invest. 1996, 56:409-14