IVD

## ANTISTREPTOLYSIN O

(ASLO 1x15)

Kat.č. 10010251

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Cat. No. 10010251

Store at (+2 to +8)°C

The liquid reagents for immunoturbidimetric quantitative determination of anti-streptolysin O in human serum.

### Principle

The human serum sample reacts upon polystyrene latex particles coated with streptolysin O. In the presence of ASO, the particles agglutinate. This induces an increase in optical density at 600 nm, which is directly proportional to the ASO titre in the sample.

### Reagents

<b>R1 Buffer</b>	(2x70 ml)
Posphate buffer pH 7	50 mmol/l
Sodium azide	0.09 g/l
<b>R2 ASLO reagent</b>	(1x15 ml)
Polystyrene particles bounds streptolysin O	
Glycine buffer pH 8.2	100 mmol/l
Sodium chloride	150 mmol/l
BSA	5 g/l
Sodium azide	0.075 g/l

### Preparation and stability of working reagents

If stored at (+2 to +8)°C, reagents are stable until expiry date, that is stated on the package. Do not freeze!

Reagents R1 and R2 are liquid ready to use. After opening, reagents are stable until expiry date at (+2 to +8)°C, if stored at appropriate temperature conditions, closed carefully and without any contamination.

Invert gently several times before use, avoiding the formation of foam.

### Samples

Serum, samples should be fresh. Do not used haemolysed or lipemic samples. The samples with particles should be centrifuged to eliminate them.

Prior the analysis, the serum can be stored 2 days at (+2 to +8)°C or 6 months at -20°C. Avoid repeated freezing and thawing.

### Calibration

For calibration use ASLO Calibrator (Cat. No. 10010252).

For photometers working with a 5 points calibration curves, prepare a range of 5 reference samples by successive 1:1 dilutions of the ASLO Calibrator by saline. If, for example, the concentration of ASO in ASLO Calibrator is 400 IU/ml, prepare the following concentrations: 25-50-100-200-400 IU/ml (dilution in ratio 1:16, 1:8, 1:4, 1:2, undiluted standard).

Plot the calibration curve.

Dilute samples with ASO concentration higher then the highest standard by saline 10x and repeat analysis. The results multiply by dilution.

### Quality control

For quality control use Bio-Rad Liquicheck™ or ILT ASLO Control (Cat. No. 10010253) or other quality control sera with ASO concentration values determined by immunoturbidimetric method.

### Procedure

Wavelength	600 nm
Cuvette	1 cm
Temperature	+37 °C
Serum/reaction mixture ratio	1/89.3

The calibrator, the controls and unknown samples must not be diluted.

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

	Sample	Standard	Reagent blank
Standard	–	0.003 ml	–
Sample	0.003 ml	–	–
Distilled water	–	–	0.003 ml
Reagent 1	0.225 ml	0.225 ml	0.225 ml
Reagent 2	0.040 ml	0.040 ml	0.040 ml

Mix R1 and R2, add sample and read turbidance A1 immediately after mixing. Allow to incubate 6 minutes at 37°C and record turbidance A2.

Applications for automatic analysers can be supplied for request.

### USED SYMBOLS

Catalogue Number

Manufacturer

See Instruction for Use

Lot Number

CE Mark - Device comply with the Directive 98/79/EC

Storage Temperature

Expiry date

In vitro Diagnostics

Content

### Calculation

$$ASO \text{ (IU/ml)} = \frac{(A_2 - A_1)_{\text{sample}} - (A_2 - A_1)_{\text{reagentblank}}}{(A_2 - A_1)_{\text{st}} - (A_2 - A_1)_{\text{reagentblank}}} \times C_{\text{cal}}$$

$C_{\text{cal}}$  ...concentration of calibrator

### Reference values

fS ASO (IU/ml) adults	< 200
fS ASO (IU/ml) children	< 100

The range of reference values is only approximate, it is recommended to all laboratories to verify the extent of the reference interval for their particular examined population.

### Performance characteristics

Limit of detection: 6 IU/ml
Low limit of quantification: 17 IU/ml
Hook effect: from 1500 IU/ml
Analytical range: 17 – 400 IU/ml

### Precision (at 37 °C)

INTRA-ASSAY n = 10	Mean (IU/ml)	SD (IU/ml)	CV (%)
Sample 1	269	8,82	3.30
Sample 2	419	8.22	1.96
INTER-ASSAY n = 10	Mean (IU/ml)	SD (IU/ml)	CV (%)
Sample 1	272	9.68	3.55
Sample 2	415	9.42	2.27

### Comparison: nephelometry

**Results:** n = 43, r = 0.997, y = 0.84 x + 17.75 IU/ml

### Interferences

There is no significant interference by bilirubin, haemoglobin or rheumatoid factor. Intralipid > 1% produces a negative interference > 10%.

### Calibration stability

On the most often automatic analysers it is > 2 weeks.

### Health protection

For in vitro diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous, but contain in low concentration (<0.1%) sodium azide, very toxic substance, dangerous for the environment.

Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.

Products containing sodium azide have to be handled with care, avoid ingestion and contact with skin and mucous membrane.

### First aid

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

### Waste disposal

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

### Reference

Difco Manual, 10th ed., Difco Laboratories, Detroit 1984, pp. 889-894.

Spaun, J. et al., International standard for anti-streptolysin-O. Bull. WHO 24 (1961) 271-279.

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