# XL ISE CLEANING

## Cleaning Solution for use with ERBA XL auto analyzers

For *in vitro* diagnostic use only.

Cat. No.	Pack Name	Packaging (Content)
XSYS0062	XL ISE CLEAN	1 x 52 ml

**INTENDED USE:** Cleaning solution is an aqueous enzyme mixture, used for preventing protein build up in the fluid lines of the ISE module in the XL Clinical Chemistry Auto-analyzer.

**WARNING:** For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagents of the kit are not classified as dangerous.

**REAGENT PREPARATION:** The reagent is ready to use.

**REAGENT COMPOSITION: Active Ingredients**

Protease (Bacterial)	0.06 µg/m <sup>3</sup>
Propylene glycol	>3%

**REAGENT STORAGE & STABILITY:** Store reagent at 15–25°C. Before use, keep the reagent at room temperature for at least 2 hours. The reagent is stable till the expiry date mentioned on the label.

**PROCEDURE:** Dispense 500 µl of cleaning solution into the XL sample cup and place on the ISE2 position of the sample tray.

The cleaning solution should be run once a day to prevent protein build up or at 8 hrs intervals if the ISE module performs more than 50 samples per day.

**PRECAUTIONS:** 1. Avoid contaminating the reagent.  
2. Keep the container tightly closed and store the reagent in the original container.

**SYMBOLS: The following symbols are used in the labelling of ERBA kits:**

REF	Catalogue Number	Manufacturer	See Instruction for Use
LOT	Lot Number	CE Mark - Device comply with the Directive 98/79/EC	Storage Temperature
Expiry Date		In Vitro Diagnostics	Content

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