EN

XL WASH

Concentrated Wash Solution

Cat. No.:	Pack Name:	Packaging (Content):
XSYS0066	XL WASH	4 x 100 ml

For *in vitro* diagnostic use only.

INTENDED USE: XL WASH is a concentrated solution used with Erba XL Automated Clinical Chemistry Analyzers.

PRINCIPLE: The performance of any Erba XL analyzer depends on the cleanliness of the microcuvetes. The use of serum or plasma or urine samples results in protein deposition. Over a period of time this may interfere with the assay resulting in erroneous results. XL WASH contains a surfactant which solubilises proteins and cleans the microcuvetes.

STABILITY AND STORAGE: The reagent is stable till the expiry date mentioned on the label when stored at room temperature (15–25°C).

Warning

Reagent 1 contains sodium hydroxide.

Hazard statement(s):

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement(s):

P302 + P352 IF ON SKIN: Wash with plenty of water and soap.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WORKING SOLUTION PREPARATION: The contents of each bottle (100 ml) of XL WASH is to be diluted to 10 litres with deionized water (10 ml XL WASH + 9900 ml de-ionized water) before using for cleaning (Cleaning solution Can) in Erba XL analyzer.

SYMBOLS:

Catalogue Number	Manufacturer	See Instruction for Use	Product name
Lot Number	CE Mark - Device comply with the Directive 98/79/EC	Storage Temperature	
Expiry Date	In Vitro Diagnostics	Content	Date of revision: 30. 7. 2015

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