# ALBUMIN

Cat. No.	Pack Name	Packaging (Content)
XSYS0001	ALB 440	R1: 10 x 44 ml

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Albumin in human serum and plasma.

## CLINICAL SIGNIFICANCE

Albumin, a major plasma protein, is synthesised in the liver from amino acids which are absorbed from the ileum. Its functions include regulation of distribution of extra-cellular fluid, transportation of various hormones, vitamins and trace metals.

### Increased levels are observed in

- Dehydration due to reduced plasma water content.
- Stasis during venipuncture which causes fluid to escape into the extravascular compartment.

### Decreased levels are observed in

- Excessive protein loss (mainly albumin) - from kidney, skin or intestine.
- Decreased synthesis due to dietary, hepatic disease or malabsorption.
- Increased catabolism in fever, untreated diabetes mellitus and hypertension.

## PRINCIPLE

Albumin binds with Bromo Cresol Green (BCG) at pH 4.2 causing a shift in absorbance of the yellow BCG dye. The blue-green colour formed is proportional to the concentration of albumin, when measured photometrically between 540-630 nm with maximum absorbance at 625 nm.

## REAGENT COMPOSITION

### R1

Bromocresol green	0.21 mmol/l
Succinate Buffer	100 mmol/l
Sodium Azide	0.5 g/l

## REAGENT PREPARATION

Reagent is liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

On board stability: min. 30 days if refrigerated (2–10°C) and not contaminated.

## SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

**Stability in serum:** 1 month at 2–8°C  
1 week at 15–25°C  
at least 3 months at -20°C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by **internal** quality control procedures

### Traceability:

This calibrator has been standardized with using the ERM-DA470.

## QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

## CALCULATION

The XL Results are calculated automatically by the instrument.

## UNIT CONVERSION

g/dl x 10 = g/l

## EXPECTED VALUES <sup>2</sup>

### Serum:

0 – 4 d	2.8 – 4.4 g/dl
4 d – 14 y	3.8 – 5.4 g/dl
14 – 18 y	3.2 – 4.5 g/dl
Adult	
20 – 60 y	3.5 – 5.2 g/dl
60 – 90 y	3.2 – 4.6 g/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.1 g/dl

**Linearity:** 6 g/dl

**Measuring range:** 0.1 – 6 g/dl

Intra-assay precision Within run (n=20)	Mean (g/dl)	SD (g/dl)	CV (%)
Sample 1	3.303	0.017	0.51
Sample 2	4.926	0.044	0.88

Inter-assay precision Run to run (n=20)	Mean (g/dl)	SD (g/dl)	CV (%)
Sample 1	3.345	0.043	1.29
Sample 2	4.612	0.043	0.93

## COMPARISON

A comparison between XL-Systems Albumin (y) and a commercially available test (x) using 40 samples gave following results:

$y = 1.000 x + 0.065$  g/dl

$r = 0.997$

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 contains < 1.2 % succinic acid.

Warning

### Hazard statement:

H319 Causes serious eye irritation.

### Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P337+P313 If eye irritation persists: Get medical advice.  
Reagent also contains less than 0.1% sodium azide - classified as very toxic and dangerous substance for the environment.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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