# ALPHA AMYLASE (SINGLE REAGENT)

Cat. No.	Pack Name	Packaging (Content)
XSYS0003	AMY 110	R1: 5 x 22 ml
XSYS0091	AMY 440 XL-1000	R1: 10 x 44 ml

EN

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of alpha-Amylase in human serum, plasma and urine.


## CLINICAL SIGNIFICANCE

$\alpha$ -Amylase is derived mainly from the salivary glands and the exocrine pancreas.  $\alpha$ -Amylase catalyses the hydrolysis of  $\alpha$ -1-4 glucosidic linkages of starch and other related polysaccharides to produce maltose and other oligosaccharides. The enzyme is a relatively small molecule which is rapidly cleared by the kidneys and excreted in the urine.

$\alpha$ -Amylase is most frequently measured in the diagnosis of acute pancreatitis when serum levels may be grossly elevated. In acute pancreatitis  $\alpha$ -amylase starts to rise approximately 4 hours after the onset of pain, reaches a peak at 24 hours and remains elevated for 3-7 days. Hyperamylasemia is also associated with other acute abdominal disorders, biliary dysfunction, salivary gland disorders, ruptured ectopic pregnancy and macroamylasemia.

## PRINCIPLE

2-Chloro-4-nitrophenol- $\beta$ -1-4 galactopyranosylmaltotrioxide (CNP-G) is a direct substrate for determination of  $\alpha$ -amylase activity, which does not require the presence of ancillary enzymes. The rate of 2-chloro-4-nitrophenol formation can be monitored at (400-420) nm and is proportional to the  $\alpha$ -amylase activity.

## REAGENT COMPOSITION

MES buffer	50 mmol/l
Calcium Chloride	3.81 mmol/l
Sodium Chloride	300 mmol/l
Potassium Thiocyanate	450 mmol/l
Sodium Azide	13.85 mmol/l
CNPG	0.91 mmol/l

## REAGENT PREPARATION

Reagent is liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2-8°C.

On board stability: min. 30 days if refrigerated (2-10°C) and not contaminated.

## SPECIMEN COLLECTION & HANDLING

Use serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

### Stability

<b>in serum/plasma:</b>	7 days	at 20-25°C
	7 days	at 4-8°C
	1 year	at -20°C
<b>in urine:</b>	2 days	at 20-25°C
	10 days	at 4-8°C
	3 weeks	at -20°C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

### Traceability:

This calibrator has been standardized to IFCC formulation.

## QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

## CALCULATION

Results are calculated automatically by the instrument.

## UNIT CONVERSION

U/l x 0.017 =  $\mu$ kat/l

## EXPECTED VALUES <sup>7</sup>

at 37°C

**Serum:** up to 80 U/l

**Urine:** up to 500 U/l

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

<b>Limit of quantification:</b>	10.8 U/l
<b>Linearity:</b>	1500 U/l
<b>Measuring range:</b>	10.8 – 1500 U/l

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	247.1	2.5	1.0
Sample 2	260.8	2.6	1.0

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	58.3	1.3	2.2
Sample 2	142.9	2.3	1.6

## COMPARISON

A comparison between XL-Systems Amylase (y) and a commercially available test (x) using 40 samples gave following results:

$y = 0.973x - 4.80$  U/l

$r = 0.989$

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 2.5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.

### Note:

Saliva and skin contain alpha-amylase therefore never pipette reagents by mouth and avoid contamination of samples and reagents. Even trace contamination can affect results.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified as dangerous but contain less than 0.1% sodium azide.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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