CALCIUM

Cat. No.	Pack Name	Packaging (Content)
XSYS0007	CA 120	R1: 10 x 12 ml
XSYS0092	CA 220 XL-1000	R1: 10 x 22 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Calcium in human serum, plasma or urine.

CLINICAL SIGNIFICANCE

Calcium has numerous function within the body, not only as a structural factor in bones and teeth, but also in normal neuromuscular function and the clotting of blood.

Hypercalcaemia may develop in patients with Paget's disease of bone and hyperparathyroidism. The cause of hypercalcaemia in malignancy is an increased bone resorption either caused by metastasis or by humoral factors produced by the tumor cell.

In Rickets, Coeliac diseases, idiopathic steatorrhea, osteomalacia, tropical sprue and following surgical resection of the small intestine, serum calcium is often moderately reduced, usually in association with low plasma protein concentration.

PRINCIPLE

Arsenazo III combines with calcium ions at pH 6.5 to form a coloured chromophore, the absorbance of which is measured at 650 nm (650-660 nm) and is proportional to calcium concentration.

Arsenazo III has a high affinity ($k^d = 1 \times 10^{-7}$) for calcium ions and shows no interference from other cations normally present in serum, plasma or urine.

REAGENT COMPOSITION

R1

Arsenazo III 0.10 mmol/l
Phosphate buffer (pH 7.8 ± 0.1) 50 mmol/l

REAGENT PREPARATION

Reagent is liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

On board stability: min. 30 days if refrigerated (2–10°C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability in serum / plasma:

7 days at 20–25°C
3 weeks at 4–8°C
8 months at -20°C

Stability in urine:

2 days at 20–25°C
4 days at 4–8°C
3 weeks at -20°C

FREEZE ONLY ONCE!

For the determination in urine dilute the sample using redistilled water in 1 + 2 ration and after adding several drops of HCl 0.1 mol/l adjust the sample pH to 3 – 4 (result x3).

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized using the SRM 956c reference material.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

mg/dl x 0.25 = mmol/l

EXPECTED VALUES ⁴

Serum:

Adult: 8.6 – 10.2 mg/dl
Child:
2 – 12 y 8.8 – 10.8 mg/dl
10 d – 24 mo 9.0 – 11.0 mg/dl
0 – 10 d 7.6 – 10.4 mg/dl

Urine:

Female < 250 mg/24 h
Male < 300 mg/24 h

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.6 mg/dl
Linearity: 16 mg/dl
Measuring range: 0.6 – 16 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	7.836	0.068	0.89
Sample 2	11.956	0.052	0.43

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	9.096	0.184	2.05
Sample 2	11.748	0.264	2.26

COMPARISON

A comparison between XL-Systems Calcium (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.979 x - 0.076 mg/dl
r = 0.999

INTERFERENCES

Following substances do not interfere:
haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 500 mg/dl.

Note:

Care must be taken to avoid calcium contamination. The use of the plastic tubes or cuvettes is strongly recommended. The user should ensure that such disposables are free for calcium contamination. If glassware is used, it should be soaked in dilute HCl deionized water and dried.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent of the kit is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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