CHLORIDE

Cat. No.	Pack Name	Packaging (Content)
XSYS0008	CL 120	R1: 10 x 12 ml

EN

CE IVD

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Chloride in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

Chloride is estimated in conjunction with the other electrolytes and used to support interpretation of the other electrolytes. The anion gap $(Na + K) - (Cl + HCO_3)$ between these electrolytes is frequently affected in disease states. Increased levels can occur in nephritis, prostatic obstruction, eclampsia and dehydration. Decreased levels can occur with impaired gastrointestinal or renal function.

PRINCIPLE

When chloride is mixed with a solution of undissociated mercuric thiocyanate, the chloride preferentially combines with mercury forming mercuric chloride. The thiocyanate that is released then combines with ferric ions present in the solution forming strongly coloured ferric thiocyanate with an absorption maxima at 480 nm.

REAGENT COMPOSITION

R1

Mercuric Thiocyanate	1.01 mmol/l
Ferric Nitrate	37.13 mmol/l
Nitric Acid	62.43 mmol/l
Mercuric Nitrate	105 µmol/l

REAGENT PREPARATION

Reagent R1 is liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

On board stability: min. 30 days if refrigerated (2–10°C) and not contaminated.

SPECIMEN COLLECTION & HANDLING

Use serum or plasma (heparin) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability:

7 days	at 20–25°C
7 days	at 4–8°C
	at least one year at -20°C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to coulometry.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

Results are calculated automatically by the instrument.

UNIT CONVERSION

mmol/l = mEq/l

EXPECTED VALUES †

Serum:

Adult	98 – 107 mmol/l
Cord	96 – 104 mmol/l
Premature	95 – 110 mmol/l
0-30 d	98 – 113 mmol/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 3.7 mmol/l

Linearity: 160 mmol/l

Measuring range: 3.7 – 160 mmol/l

Intra-assay precision Within run (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	116.1	1.25	1.08
Sample 2	108.3	1.64	1.52

Inter-assay precision Run to run (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	117	1.47	2.75
Sample 2	94.6	1.43	1.52

COMPARISON

A comparison between XL-Systems Chloride (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 1.000 x + 1.0 \text{ mmol/l}$$

$$r = 0.927$$

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 500 mg/dl.

Note:

Grossly haemolysed, lipaemic and jaundiced samples give falsely elevated results and should not be used. High levels of immunoglobulins, as found in case of multiple myeloma will interfere due to the development of turbidity.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 contains < 6% methanol.

Warning

Hazard statement:

H302 Harmful if swallowed.

H312 Harmful in contact with skin.

H332 Harmful if inhaled.

H371 May cause damage to organs.

Precautionary statement:

P260 Do not breathe vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection.

P301 + P312 IF SWALLOWED: Call a doctor if you feel unwell.

P302 + P352 IF ON SKIN: Wash with plenty of water and soap.

WASTE MANAGEMENT

Please refer to local legal requirements.

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