# CREATININE

Cat. No.	Pack Name	Packaging (Content)
XSYS0024	CREA 275	R1: 5 x 44 ml, R2: 5 x 11 ml
XSYS0076	CREA 564 XL-1000	R1: 6 x 72 ml, R2: 6 x 22 ml

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Creatinine in human serum, plasma and urine.

## CLINICAL SIGNIFICANCE

Creatinine is a waste product formed in muscle from the high energy storage compound, creatine phosphate. The amount of creatinine produced is fairly constant (unlike Urea) and is primarily a function of muscle mass. It is not greatly affected by diet, age, sex or exercise. Creatinine is removed from plasma by glomerular filtration and then excreted in urine without any appreciable resorption by the tubules.

Creatinine is used to assess renal function, however, serum creatinine levels do not start to rise until renal function has decreased by at least 50%.

## PRINCIPLE

Creatinine reacts with alkaline picrate to produce a reddish colour (Jaffe reaction). This is a non-specific reaction and is given by many other substances. Specificity of the assay has been improved by the introduction of a kinetic method<sup>1</sup>, however, the cephalosporin antibiotics are still major interferents.

## REAGENT COMPOSITION

<b>R1</b>	
Sodium Hydroxide	240 mmol/l
<b>R2</b>	
Picric Acid	26 mmol/l

## REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

On board stability: min. 12 days if refrigerated (2–10°C) and not contaminated and if chimneys are used.

Absorption of atmospheric CO<sub>2</sub> in open reagent bottle R1 leads to the reduced stability of reagents. To minimize the rate of absorption of CO<sub>2</sub> the use of “chimneys” provided with the kit is recommended. Use of chimneys has shown improvement in calibration stability, provided the storage and use of reagents is as per instructions.

Chimneys provided with the kit are to be placed directly into the reagent bottle immediately after opening the cap.

White chimney is to be used for reagent R1 and black chimney for reagent R2. Color differentiation of chimneys is just to ensure chimneys are not exchanged between R1 & R2.

Chimneys can be re-used again for the bottle of the same kit.

Before using the chimney in a new bottle, it is recommended to wash chimneys well with distilled water and dry them thoroughly.

## SPECIMEN COLLECTION & HANDLING

Use unheamolytic serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

### Stability

<b>in serum /plasma:</b>	7 days	at 4–25°C
	at least 3 months	at -20°C
<b>in urine:</b>	2 days	at 20–25°C
	6 days	at 4–8°C
	6 months	at -20°C

For the determination in urine use 24 hours specimen. It is important to exactly measure the volume of collected urine. Dilute urine samples in 1+19 ratio with distilled water and multiply results by 20.

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

### Traceability:

This calibrator has been standardized to ID-MS.

## QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

## CALCULATION

Results are calculated automatically by the instrument.

## UNIT CONVERSION

mg/dl x 88.4 = μmol/l

## EXPECTED VALUES <sup>3</sup>

	Serum	Urine
Males:	0.7 - 1.3 mg/dl	14 - 26 mg/kg/day
Females:	0.6 - 1.1 mg/dl	11 - 20 mg/kg/day
Newborn:	0.3 - 1.0 mg/dl	
Infant:	0.2 - 0.4 mg/dl	8 - 20 mg/kg/day
Child:	0.3 - 0.7 mg/dl	8 - 22 mg/kg/day
Adolescent:	0.5 - 1.0 mg/dl	8 - 30 mg/kg/day

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.08 mg/dl

**Linearity:** 18 mg/dl

**Measuring range** 0.08 - 18 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	2.97	0.043	1.45
Sample 2	4.49	0.052	1.16

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	1.25	0.021	1.71
Sample 2	3.31	0.031	0.95

## COMPARISON

A comparison between XL-Systems Creatinine (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 0.988x + 0.029 \text{ mg/dl}$$

$$r = 0.997$$

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 15 mg/dl, triglycerides up to 1000 mg/dl.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 contains 1.0 % sodium hydroxide.

Warning

### Hazard statement:

H315 Causes skin irritation

H319 Causes serious eye irritation

### Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Fabiny DL and Ertinghausen G. Clin Chem, 1971 : 71 : 391.
2. Kroll MH and Elin R.J. Clin Chem, 1983 : 29 : 2044.
3. Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
4. Young DS, Effects of Drugs on Clinical Laboratory Tests. Third Edition 1990; 3 : 122-32.
5. Henry RJ (Ed), Clinical Chemistry : Principles and Technics (2nd Ed), Clinical Chemistry : Principles and Technics (2nd Ed), Harper and Row, 1974; 548-551.
6. Wachtel et al, Creation and Verification of Reference Intervals. Laboratory Medicine 1995; 26-593-7.
7. National Committee for Clinical Laboratory Standards. User Evaluation of Precision Performance of Clinical Chemistry Devices. NCCLS; 1984, NCCLS Publication EP5-T.

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