UIBC

Cat. No.	Pack Name	Packaging (Content)
XSYS0050	UIBC 125	R1: 4 x 25 ml, R2: 2 x 12.5 ml, R3 STD: 2 x 2 ml
XSYS0090	UIBC 220 XL-1000	R1: 4 x 44 ml, R2: 2 x 22 ml, R3 STD: 2 x 2 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of unsaturated iron-binding capacity in serum and plasma. Photometric test using chromogen ferrozine.

CLINICAL SIGNIFICANCE

Iron exists in serum complexed with transferrin, a transport protein. Most early procedures for iron determination involved dissociation of the iron from the iron-protein complex, precipitation of the proteins, and then measurement of the iron content of the protein free filtrate.

Many chromogens have been used in the determination including thiocyanate o-phenantroline, bathophenantroline and TPTZ. In 1971 Presijn et al.¹ presented a method using the chromogen ferrozine, described by Stookey.² This method did not require protein precipitation and was more sensitive than previous methods.

The present procedure is a modification of the Presijn method.

In most cases, both serum iron and TIBC values are necessary for greatest diagnostic significance. Low serum iron values are seen in chronic blood loss, insufficient intake or absorption of iron and increased demand on the body stores (e.g. pregnancy). Elevated serum iron values are seen in increased red cell destruction, decreased red cell synthesis, increased iron take, or increased iron stores release.

Increase in the TIBC may be due to increased production of apotransferrin (e.g. chronic iron deficiency) or an increased release of ferritin, as in hepatocellular necrosis.

Decreases in the TIBC can occur with cirrhosis and hemochromatosis due to a deficiency in ferritin, or in nephrosis due to a loss of apotransferrin.

PRINCIPLE

Total Iron-Binding Capacity (TIBC): A known amount of ferrous ions are added to serum at an alkaline pH. The ferrous ions bind with transferrin at unsaturated iron-binding sites. The additional unbound ferrous ions are measured using the ferrozine reaction. The difference between the amount of ferrous ions added and the unbound ions measured is the unsaturated iron-binding capacity (UIBC). The TIBC is equal to the serum iron concentration plus the UIBC.

REAGENT COMPOSITION

R1

Tris buffer (pH 8.45)	220 mmol/l
Ferrous ammonium sulfate	12.1 µmol/l
Hydroxylamine hydrochloride	100 mmol/l

R2

Hydroxylamine hydrochloride	220 mmol/l
Ferrozine	≥ 3.0 mmol/l

R3 STD

Iron standard	500 µg/dl (89.5 µmol/l)
---------------	-------------------------

REAGENT PREPARATION

All reagents are ready to use.

STABILITY AND STORAGE

Unopened kit components up to stated expiration date at 2–8°C.

Reagents are light-sensitive. Do not let bottles remain open. Keep containers tightly closed.

SPECIMEN COLLECTION AND HANDLING

Use serum or lithium heparinized plasma not haemolyzed. Separate serum from clot within one hour after collection. Separate plasma from cells within one hour after collection.

It is recommended to follow NCCLS procedures (or similar standardized conditions)

Stability: 7 days at 4°C
4 days at 20–25°C

Centrifuge samples containing precipitates before performing the assay.

CALIBRATION

For calibration it is recommended the standard included in the set.

QUALITY CONTROL

For quality control Erba Norm, Cat. No. BLT00080, Erba Path, Cat. No. BLT00081 are recommended.

CALCULATION

Results are calculated automatically by the instrument.

UIBC Calculation

Iron Level + UIBC = TIBC (µg/dl)

UNIT CONVERSION

µg/dl x 0.179 = 1 µmol/l

EXPECTED VALUES ³

UIBC: 110 – 370 µg/dl	UIBC: 19.7 – 66.2 µmol/l
TIBC: 228 – 428 µg/dl	TIBC: 40.8 – 76.6 µmol/l
(TIBC = Iron + UIBC)	

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	12.2 µg/dl
Linearity:	830 µg/dl
Measuring range:	12.2 – 830 µg/dl (2.18 – 148.6 µmol/l)

PRECISION

Intra-assay precision Within run (n=20)	Mean (µg/dl)	SD (µg/dl)	CV (%)
Sample 1	366	16.8	4.58
Sample 2	487	11.5	2.35

Inter-assay precision Run to run (n=20)	Mean (µg/dl)	SD (µg/dl)	CV (%)
Sample 1	454	16.8	3.69
Sample 2	234	10.2	4.37

COMPARISON

A comparison between XL-Systems UIBC Ferrozine (y) and a commercially available test (x) using 40 samples gave following results.

r = 0.990 y = 1.046 x – 2.626 µg/dl

INTERFERENCES

Following substances do not interfere:
haemoglobin up to 100 mg/dl, bilirubin up to 20 mg/dl, triglycerides up to 1250 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent R1 contains <0.7 % hydroxylamine hydrochloride and may produce an allergic reaction.

Reagent R2 contains less than 0.1 % sodium azide - classified as very toxic and dangerous substance for the environment.

Reagent R2 and R3 STD contain 1.5 % hydroxylamine hydrochloride.

R2 and R3 STD:

Warning

Hazard statement:

H317 May cause an allergic skin reaction.

H351 Suspected of causing cancer.

Precautionary statement:

P202 Do not handle until all safety precautions have been read and understood.

P261 Avoid breathing vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P308 + P313 IF exposed or concerned: Get medical advice.

P333 + P313 If skin irritation or rash occurs: Get medical advice.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	UIBC	UIBC	UIBC	UIBC	UIBC	UIBC
Test Code	58	58	58	58	58	58
Report Name	UIBC	UIBC	UIBC	UIBC	UIBC	UIBC
Unit	µg/dl	µg/dl	µg/dl	µg/dl	µg/dl	µg/dl
Decimal Places	1	1	1	1	1	1
Wavelength-Primary	578	578	570	570	570	578
Wavelength-Secondary	700	700	700	700	700	700
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	Linear	Linear	Linear	Linear	Linear	Linear
M1 Start	16	16	12	22	10	16
M1 End	16	16	12	22	10	16
M2 Start	34	36	51	63	29	34
M2 End	34	36	51	63	31	34
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	2.5	2.5	2.5	2.5	2.5	2.5
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	12.2	12.2	12.2	12.2	12.2	12.2
Technical Maximum	830	830	830	830	830	830
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	0	0	0	0	0	0
Reagent Abs Max	0.2	0.2	0.2	0.2	0.2	0.2
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	UIBC R1	UIBC R1	UIBC R1	UIBC R1	UIBC R1	UIBC R1
Reagent R2	UIBC R2	UIBC R2	UIBC R2	UIBC R2	UIBC R2	UIBC R2
Reagent R3	NA	NA	NA	NA	NA	NA

Test Volumes						
Test	UIBC	UIBC	UIBC	UIBC	UIBC	UIBC
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes						
Normal	15	15	15	15	15	15
Dilution Ratio	1	1	1	1	1	1
Increase	30	30	30	30	30	30
Dilution Ratio	1	1	1	1	1	1
Decrease	5	5	5	5	5	5
Dilution Ratio	1	1	1	1	1	1
Standard volume	15	15	15	15	15	15
Reagent Volumes and Stirrer speed						
RGT-1 Volume	200	200	200	200	160	200
R1 Stirrer Speed	High	High	NA	High	High	High
RGT-2 Volume	50	50	50	50	40	50
R2 Stirrer Speed	High	High	NA	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA

Reference Ranges						
Test	UIBC	UIBC	UIBC	UIBC	UIBC	UIBC
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	110	110	110	110	110	110
Normal-Upper Limit	370	370	370	370	370	370
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	110	110	110	110	110	110
Normal-Upper Limit	370	370	370	370	370	370
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA

Revision Number						
Revision	<A-100- UIBC-1 20.08.2013>	<A-200- UIBC-1 20.08.2013>	<A-300/600- UIBC-1 20.08.2013>	<A-640- UIBC-1 20.08.2013>	<A-1000- UIBC-1 20.08.2013>	<A-180- UIBC-1 12.12.2013>

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1. Goodwin J, Murphy D, Guillemette M: Clin Chem 1966;12:47
2. Henry RJ, Cannon DC, Winkelman W: Clinical Chemistry Principles and Techniques. Hagerstown, MD, Harper & Row, Inc, 1974; 684.
3. Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, Pa:WB Saunders Co ; 1995; 376.
4. Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements, from two different analytical methods. J Clin Chem Biochem 1983; 21:709-720.
5. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Biochem 1988; 26:783-790.
6. Persijn, J.P., et al. Clin. Acta 35.91 (1971).
7. Stookey, L.L., Anal. Chem 42:779 (1970).
8. Weissman, N., Pileggi, V.J., in Clinical Chemistry: Principles and Technics, 2nd Ed., R.J.Henry et al, editors, Hagerstown (MD), Harper & Row, pp. 692-693 (1974).
9. Young, D.S. et al, Clin Chem. 21: 1D (1975)
10. Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods, Philadelphia, W.B. Saunders, P.1434 (1984).

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