IVD

Činidlo pro test PYR

CZ

Kat. č.: 10003379

Pro mikrobiologii

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Cat. No.: 10003379

For microbiology

Reagent for PYR test is the colour forming reagent for performing of pyrrolidonylarylamidase test (PYRAtest strips and/or test PYR in MIKROLATEST® kits).

Preparation for use:

Pour the content of the bottles labelled Reagent for PYR I and Reagent for PYR II into the glass bottle labelled Reagent for PYR and mix well. The dropping bottle is intended for adding the Reagent for PYR into wells of microplate. A little pipette included in the PYRAtest kit is intended for adding of the reagent on the strip zone.

Reagent for PYR test contains:

Reagent for PYR I	9 ml
Reagent for PYR II	9 ml
Dropping bottle of Reagent for PYR	1 pc
Information leaflet	

Storage, expiration:

Reagent for PYR test and mixed Reagent for PYR should be stored at the temperature of (+2 to +8) °C until the expiry date indicated on the package.

Note: in the solution of Reagent for PYR I and in the mixed working solution of Reagent for PYR may after prolonged storage to form a suspension which is not harmful.

Warning:

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent for PYR I (35%) and Reagent for PYR (17%) contain harmful 2-methoxyethanol, which is toxic for fertility (category 1B). If swallowed it is necessary to drink about 0.5 l of lukewarm water. On skin or eye contact rinse the place with jet of tap water.

If necessary, consult a physician.

Danger

Hazard statement:

H302 Harmful if swallowed.

H312 Harmful in contact with skin.

H332 Harmful if inhaled.

H360FD May damage fertility. May damage the unborn child.

Precautionary statement:

P201 Obtain special instructions before use.

P280 Wear protective gloves/protective clothing/eye protection.

P308+P313 IF exposed or concerned: Get medical advice.

Waste disposal:

Contingent rest of the reagent should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials.

Put packaging paper waste and containers to recycling.

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