## Opti-sol XL 750

Cat. No.:	Pack Name:	Packaging (Content):
REG00054	Opti-sol XL 750	6 x 100 ml

For *in vitro* diagnostic use only.

The part of set is chip RFID.

**Intended use:** Opti-sol XL 750 is a concentrated solution for use with fully automated urine analyser LAURA XL.

**Principle:** The performance of analyser LAURA XL depends on stabilization and cleanness of purity optical and flow parts. Opti-sol XL 750 solution contains additives which stabilize optical and fluidic parts and provides clean environment in the system.

**Stability and storage:** The reagent is stable till the expiry date mentioned on the label, when stored at room temperature (15–25 °C).

The solution contains sodium hypochlorite.

Danger

**Hazard statement(s):**

H314 Causes severe skin burns and eye damage. H335 May cause respiratory irritation.  
H400 Very toxic to aquatic life. EUH031 Contact with acids liberates toxic gas.

**Precautionary statement(s):**

P273 Avoid release to the environment.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing.  
Rinse skin with water/shower.  
P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

**Working solution preparation:** The content of the bottle of Opti-sol XL 750 (100 ml) dilute in 5 litres of deionized water (100 ml of Opti-sol XL 750 + 4900 ml of deionized water) in to provided can and attach to LAURA XL instrument. Working solution is not classified as dangerous.

**Waste disposal:** Residues of chemicals and empty containers of chemicals should be disposed of in accordance with applicable laws.

Used symbols:

Catalogue Number

Manufacturer

See Instruction for Use

Opti-sol XL 750

Product name

Lot Number

CE Mark - Device comply with the Directive 98/79/EC

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