Opti-sol XL 1500

Cat. No.:	Pack Name:	Packaging (Content):
REG00055	Opti-sol XL 1500	1 x 200 ml

For *in vitro* diagnostic use only.

The part of set is chip RFID.

Indented use: Opti-sol XL 1500 is a concentrated solution for use with fully automated urine analyser LAURA XL.

Principle: The performance of analyser LAURA XL depends on stabilization and cleanest purity of optical and flow parts.

Opti-sol XL 1500 solution contains additives which stabilize optical and fluidic parts and provides clean environment in the system.

Stability and storage: The reagent is stable till the expiry date mentioned on the label, when stored at room temperature (15–25 °C).

The solution contains sodium hypochlorite.

Danger

Hazard statement(s):

H314 Causes severe skin burns and eye damage
H400 Very toxic to aquatic life.

H335 May cause respiratory irritation.
EUH031 Contact with acids liberates toxic gas.

Precautionary statement(s):

P273 Avoid release to the environment.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Working solution preparation: The content of the bottle of Opti-sol XL 1500 (200 ml) dilute in 10 litres of deionized water (200 ml of Opti-sol XL 1500 + 9800 ml of deionized water) in to provided can and attach to LAURA XL instrument. Working solution is not classified as dangerous.

Waste disposal: Residues of chemicals and empty containers of chemicals should be disposed of in accordance with applicable laws.

Used symbols: Catalogue Number Manufacturer See Instruction for Use Product name
 Lot Number CE Mark - Device comply with the Directive 98/79/EC Storage Temperature Date of revision: 2.11. 2016
 Expiry Date In Vitro Diagnostics Content

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