IVD

OPTOCHIN

Kat. č.: 50001522

Pro mikrobiologii

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Kat. č.: 50001522

For microbiology

Diagnostic discs OPTOCHIN are intended for simple routine presumptive detection of Streptococcus pneumoniae.

Principle:

Test is based on the susceptibility of Streptococcus pneumoniae to optochin (ethylhydrocuprein hydrochlorid). Streptococcus pneumoniae shows after incubation defined zone of inhibition around the impregnated disc, whilst other viridans-group streptococci are resistant to optochin or show a very small zone of inhibition. Disc OPTOCHIN may also be applied on primary culture plate in order to provide rapid indication of presence of large numbers of Streptococcus pneumoniae in sample.

A package of the OPTOCHIN contains:

- 100 diagnostic discs for 100 determinations
- working leaflet

Storage, expiration:

Discs OPTOCHIN should be stored in dark in a refrigerator at (+2 to +8)°C until the expire date indicated on the package.

Required material:

- Blood agar plates
- common laboratory equipment (loops, sterile forceps, burner)

Recommended procedure:

- Prepare suspension corresponding to 0,5 McFarland value in saline (0,85% NaCl) from 24 hours old culture. Streak with a swab on a blood agar plate to obtain almost confluent growth.
- Alternative procedure: Pour the prepared suspension on the blood agar plate entirely, suck-off the excess of suspension by pipette and leave the agar plate to dry up for cca 20 min.

Note: Simultaneously providing of positive control with known strain is recommended.

- Aseptically place disc OPTOCHIN on the surface of the medium
- Incubate in an atmosphere containing (5–10)% CO₂ at 37°C for (18–24) h.

Caution:

- For professional use only.
- After use all discs must be autoclaved or incinerated.

Evaluation:

After incubation read the diameters of inhibition zones:

Inhibition zone (mm)	Result
> 14	Streptococcus pneumoniae
0	viridans-group streptococci other than Streptococcus pneumoniae
< 14	Streptococcus pneumoniae must be confirmed with bile solubility test

Standard control strains:

Quality control of diagnostic discs OPTOCHIN may be provided with standard control strains:

Strain	Control
Streptococcus pneumoniae CCM 4501	positive
Streptococcus constellatus CCM 4043	negative

These strains are distributed by:

CCM – Czech Collection of Microorganisms, Masaryk University, Faculty of Science, Kamenice 5, building A25, 625 00 Brno, CZ.
Tel. +420 549 491 430, Fax +420 549 498 289, <http://www.sci.muni.cz/ccm>, e-mail: ccm@sci.muni.cz
Strains are supplied in freeze-dried form or on gelatin discs.

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