Nr kat.: 10020284

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Cat. No.: 10020284

For microbiology

The kit is designed to test antimicrobial susceptibility of bacteria isolated from urine tract, especially family *Enterobacteriaceae*, on the basis of bacterial growth at breakpoint concentrations. The kit contains 40 detection strips.

Breakpoint concentrations and interpretation of sensitivities were determined on the basis of the European standard EUCAST ([www.eucast.org](http://www.eucast.org)) dated from 5.1.2011. The test is based on the rehydration of antibiotics in the wells with Suspension medium for SENSILATEST and addition of bacterial suspension. The results are read visually after 18-24 hours of incubation. We recommend to use other kits (SENSILAtest NEFERM, SENSILAtest STAPHY, SENSILAtest G+) for other bacteria group than *Enterobacteriaceae*.

**The kit contains:**

- 10 plates for examination
- A lid
- 10 pc of PE bags
- Record sheet
- Storage bag
- 1 pc of Suspension medium for SENSILATEST in dehydrated form – Attention: It is necessary to prepare it before testing!

To avoid preparation of suspension media from dehydrated form, it is possible to order ready to use media. Ready to use media are available as item 10020286 - Suspension medium for SENSILATEST. Each kit contains 20 glass tubes. One tube contains 4.4 ml of suspension media, which allows one examination on any type of SENSILATEST strip.

**Storage and expiration of the kit:**

It is recommended to store the kit at (+15 to +25)°C. The date of expiration is indicated on each package. After the aluminium package is opened, it is possible to insert any unused strips into a storage bag with desiccant at the room temperature for a period of maximum 8 days. Mark the storage bag with a type of stored strips (e.g. G-I, G-II etc.) to avoid mistake in their further usage. **Don't leave opened strips unprotected!!! Exposure to open air leads to antibiotic activity failure!!!** In case that laboratory is not able to keep kits at recommended room temperature (not exceeding 25 °C), it is better to keep kits in refrigerator (+2 to +8) °C. Then it is necessary to leave plate at room temperature for 2 hours before the use to avoid water condensation.

**Storage and expiration of Suspension medium for SENSILATEST:**

It can be stored in powder form at (+15 to +25)°C till the date indicated on the package. The expiration of the Suspension medium for SENSILATEST is 6 months at (+15 to +25)°C after it was prepared into a liquid form. Do not store the prepared medium at lower temperatures as it can lead to formation of gel in the media. If such a situation occurs, adjust the medium to room temperature and shake it thoroughly before the use.

**Material required to perform a test, not included in the kit:**

- Sterile physiological solution
- Sterile tubes
- An stepper pipette for dosage of 100 µl
- A pipette for dosage of 20 µl
- Densitometer (e.g. Densi-La-Meter II, Erba Lachema)
- Incubator 37 °C

**Material required to prepare Suspension medium for SENSILATEST, not included in the kit:**

- A screw top bottle
- A mixing stick or a stirrer
- Autoclave

**Caution:** *The kit is for professional use only! Respect the rules for work with infectious material!*

### Instructions for Use

**Caution:** *Prepare the Suspension medium for SENSILATEST before you start testing!*

To avoid preparation of suspension media from dehydrated form, it is possible to order ready to use media. Ready to use media are available as item 10020286 - Suspension medium for SENSILATEST. Each kit contains 20 glass tubes. One tube contains 4.4 ml of suspension media, which allows one examination on any type of SENSILATEST strip.

**Preparation of Suspension medium for SENSILATEST:**

One bag is sufficient to prepare 0.5 L of Suspension medium for SENSILATEST. Pour the bag content in to a screw top bottle. Add 0.5 L of preheated distilled water and stir it. Autoclave the prepared suspension for 30 minutes at 121 °C. Cool it to room temperature and stir it again under sterile conditions. A small precipitate at the bottom does not influence the function of the medium. Check the pH of the medium before the use. pH of the medium must be 7,5 (± 0,1). Adjust pH if necessary.

**Instructions to use the kit:**

- 1) Prepare a tube with 2 - 5 ml of physiological solution of pH 5,5 – 6,5.
- 2) Prepare a tube with Suspension medium for SENSILATEST optimally preheated to 37 °C.
- 3) Remove few colonies from 18 – 24 hour culture on blood agar and prepare a bacterial suspension of density of 0.5 on McF scale in physiological solution.  
**Overestimation of suspension density is smaller error than its underestimation.**
- 4) Inoculate 20 µl of bacterial suspension into 4.4 ml of Suspension medium for SENSILATEST and homogenise well.

**Inoculation:**

Remove a plate from aluminium bag. Remove required number of strips for examination. Do not remove the strips from the package earlier than 30 minutes before you start the work. Insert any unused strips into a ZIP log bag. Mark the storage bag with a type of stored strips (e.g. G-I, G-II etc.) to avoid mistake in their further usage. Transfer a desiccant bag from the original package into a ZIP log bag and close it for later use. Transfer a desiccant bag from the original package into a ZIP log bag and close it for later use. Use the strips stored in this way within 8 days! **Don't leave opened strips unprotected!!! Exposure to open air leads to antibiotic activity failure!!!**

Record numbers of the examined strains on the corresponding strips. Mark the frame with a type of strip (e.g. G-I, G-II etc.) to avoid mistake in reading results after incubation.

Inoculate 100 µl of suspension in Suspension medium for SENSILATEST into each well of the strip with the help of automatic stepper pipette.

**Incubation:**

Insert the inoculated plate into a PE bag. Fold the open end of the bag under the plate to prevent evaporation during the incubation. Incubate the plate at 37 °C for 18 – 24 hours.

**Evaluation:**

Remove the plate from the PE bag. To read the growth in the microwells, choose a way which is the most convenient for you:

- 1) Use the lid with imprinted antibiotics sequence.
- 2) Read against a grey background.
- 3) Read against natural or artificial dispersed light – you can detect even limited growth in a form of grains or light turbidity.

# MIKROLATEST®

**Please read with attention:**

You must see a growth in a control well! If the growth is not present, the test **MUST NOT** be evaluated! However, the growth control is not used to compare the presence of the growth in individual wells with antibiotics. Every change from a total transparency is evaluated as growth, including small grains of growth in a well, light turbidity or partial growth covering only a part of a well. Beware to differentiate grains of growth from media bubbles. Record the results.

H	G	F	E	D	C	B	A
Growth control	Ampicillin	Amoxicillin / clavulanate	Piperacillin / tazobactam	Cefotaxime		Ceftazidime	
GC	AMP 8	AMC 8/4	PIT 16/4	CTX 1	CTX 2	CAZ 1	CAZ 4
Ceftazidime / clavulanate		Gentamicin		Amikacin		Norfloxacin	
CZC 1/4	CZC 4/4	GEN 2	GEN 4	AMK 8	AMK 16	NOR 0,5	NOR 1
Ciprofloxacin		Trimethoprim / sulfamethoxazol		Ertapenem		Fosfomycin	Nitrofurantoin
CIP 0,5	CIP 1	T/S 2/38	T/S 4/76	ERT 0,5	ERT 1	FOS 32	NFT 64

**Interpretation:**

Growth control	GC	+										Results: +/-	Evaluation (S, I, R)
		S	I	R	S	I	R	S	I	R	S		
Amikacin	AMK	-	-	S	+	-	I	+	+	R			
Ampicillin	AMP	-	-	S				+		R	x		
Amoxicillin/clavulanate	AMC	-	-	S				+		R	x		
Ceftazidime	CAZ	-	-	S	+	-	I	+	+	R			
Ceftazidime / clavulanate	CZC	-	-	S	+	-	I	+	+	R			
Cefotaxime	CTX	-	-	S	+	-	I	+	+	R			
Ciprofloxacin	CIP	-	-	S	+	-	I	+	+	R			
Ertapenem	ERT	-	-	S	+	-	I	+	+	R			
Fosfomycin	FOS			S				+		R	x		
Gentamicin	GEN	-	-	S	+	-	I	+	+	R			
Nitrofurantoin	NFT	-	-	S				+		R	x		
Norfloxacin	NOR	-	-	S	+	-	I	+	+	R			
Piperacillin / tazobactam	PIT	-	-	S				+		R	x		
Trimethoprim / sulfamethoxazol	T/S	-	-	S	+	-	I	+	+	R			

If there is a growth at a higher concentration of an antibiotic and no growth at a lower concentrations of the same antibiotic, do not evaluate susceptibility to this antibiotic. This combination shows an error in the procedure!

Resistance to **CTX**, **CAZ** indicates possibility of **ESBL** production. It is necessary to perform disk diffusion test with disks cefpodoxim, cefpodoxime/clavulanate event. amoxicillin, amoxicillin/clavulanate and compare size of inhibition zones.

It is necessary to compare results of sensitivity to **NOR** and **CIP**. Result **NOR sensitive** along with **CIP resistant is very improbable (false)** and should be considered as an **error** in test procedure.

**Resistance to ERT** indicates **resistance to carbapenems** and it is necessary to test IMI and MER by disk diffusion test or by gradient MIC test(e.g. E-test).

Result **AMP sensitive** along with **AMC resistant** is improbable.

Depending on national or laboratory standards, other criteria described in EUCAST Expert Rules in Antimicrobial Susceptibility Testing can be used.

**Quality control:**

Control strains can be used for internal testing of functionality of the antibiotics in the laboratory. Please refer to the acceptable limits for quality control strains according to EUCAST.

CCM 3954 <i>Escherichia coli</i> (ATCC 25922)							
GC +	AMP - 8	AMC - 8/4	PIT - 16/4	CTX - 1	CTX - 2	CAZ - 1	CAZ - 4
CZC <sup>1</sup> - 1/4	CZC <sup>1</sup> - 4/4	GEN - 2	GEN - 4	AMK - 8	AMK - 16	NOR - 0,5	NOR - 1
CIP - 0,5	CIP - 1	T/S - 2/38	T/S - 4/76	ERT - 0,5	ERT - 1	FOS - 32	NFT - 64

1 no result stated in EUCAST, result stated by internal quality control ERBA Lachema s.r.o.

+ growth

- no growth

**Other recommended control strains:**

CCM 3955 *Pseudomonas aeruginosa* (ATCC 27853)

CCM 4223 *Staphylococcus aureus* (ATCC 29213)

ATCC – American Type Culture Collection

CCM – Czech Collection of Microorganisms

**Health protection:**

Components of the kit do not contain any dangerous substances.

**Disposal of the used material:**

Insert the used plate into the vessel intended for the infectious material and autoclave or destroy it by incineration.

Put paper packaging waste to recycling.

Date of last revision: 24.1.2014